

PRODUCT CATALOG 2020

bip MED[®]
hospital equipments solutions
PT. BUMI INDAH PUTRA



Management System
ISO 9001:2015
www.tuv.com
ID 9108627873



SERTIFIKAT
CPAKB
NO. FK.01.04/VI/953/2018





PT. BUMI INDAH PUTRA



PRODUCT CATEGORY



Product Name	Type	Page
COVID 19 ISOLATION CHAMBER		9
Disinfectant Chamber	BDC-102-DC	
Swap Chamber	BDC-103-SC	
Isolation Chamber	BAS-702-IC	
Isolation Chamber	BAS-502-IS	
AMBULANCE		10
Ambulance	Deluxe	
Ambulance	Standart	
Ambulance	Ekonomi	
Ambulance	Jenazah	
AMBULANCE STRETCHER		15
Ambulance Stretcher Standard (Hollow Alumunium)	BAS-302-S	
Ambulance Stretcher Multipurpose	BAS-403- MS	
Ambulance Stretcher folding standard	BAS- 401-SS	
Ambulance Stretcher (Tandu S/S) / Folding Stretcher	BAS-201-S	
BABY CARE		17
Baby Box P/C	BBB- 401	
Baby Box S/S	BBB-201	
Baby Cot P/C / Baby Basket	BBC-401-S	
Baby Cot S/S / Baby Basket	BBC-201-R	
Baby Tafel (Resussitasi Table) P/C	BTA-101	
Baby Tafel (Resussitasi Table) S/S	BTA-201	
Baby Tafel (Resussitasi Table) + Stand Lamp P/C	BTA-101-L	
Baby Tafel (Resussitasi Table) + Stand Lamp S/S	BTA - 401 - LS	
Baby Tripod (Timbangan Posyandu)	BBT-301	
Baby Measurement / Alat ukur bayi alumunium	BBT-302-A	
Baby Measurement / Alat ukur bayi Kayu	BBT-302-K	
BED SCREEN		21
Bed Screen Single P.C	BSS-101-SR	
Bed Screen Single S/S	BSS-201-SR	
Bed Screen Double P/C	BSS-102-DR	
Bed Screen Double S/S	BSS-202-DR	
Bed Screen Triple P.C	BSS-103-TR	
Bed Screen Triple S/S	BSS-203-TR	

Product Name	Type	Page
BED SIDE CABINET		22
Bed Side Cabinet (Mimbar)	BSC-101	
Bed Side Cabinet (Food Tray Mimbar)	BSC-102-FT	
Bed Side Cabinet ABS	BSC-101-A	
Bed Side Cabinet	BSC-104	
Bed Side Cabinet	BSC-105-TO	
BLUE LIGHT PHOTO THERAPY		23
Blue Light / Photo Therapy 4 Lamp	BPT-301-4L	
Blue Light / Photo Therapy 5 Lamp	BPT-301-5L	
Blue Light / Photo Therapy 6 Lamp	BPT-301-6L	
BIOLOGICAL SAFETY CABINET		24
Biological Safety Cabinet	BSF-101-C	
Biological Safety Cabinet	BSF-101	
BRANGCART PATIENT / TRANSFER PATIENT		25
Brangcart Patient P/C	BBP-101	
Brangcart Patient S/S	BBP - 401 - S	
Brangcart Standard P/C	BBP-102-S	
Brangcart Standard S/S	BBP-202-S	
Brangcart Complete P/C	BBP-103-C	
Brangcart Complete S/S	BBP-403 - CS	
Brangcart Complete Hi-Lo P/C	BBP- 404 - HLR	
Brangcart Super Deluxe	BBP-111 - DA	
Transfer Stretcher Trolley	BBP-110-TST	
CHAIR		28
Phlebotomi Chair P/C (Kursi Transfusi Darah)	BPC-101	
Phlebotomi Chair S/S (Kursi Transfusi Darah)	BPC-201	
Phlebotomi Chair Electric (Kursi Transfusi Darah)	BPC-201-E	
Hemodialysis Chair	BKR-101-HM	
Wheel Chair	BKR-202-R	
Wheel Chair	BKR-203-PP	
EXAMINATION TABLE		30
Examine Table P/C (Hollow) + Foot Step	BET-401-PS	
Examine Table (Pipa S/S) + Foot Step	BET-401-SS	
Examine Table Advance + Foot Step	BET-102-A	
Examine Table Electric P/C + Foot Step	BET-101-E	
Examine Table Electric S/S + Foot Step	BET-201-E	

Product Name	Type	Page
EXAMINE LAMP / STAND LAMP		31
Examine Lamp Halogen	BEL-102-H	
Examine Lamp LED	BEL-102-LED	
Examine Lamp LED	BEL-104-LD	
FOOT STEP		32
Foot Step P/C	BFS-101	
Foot Step S/S	BFS-201	
GYNAECOLOG		33
Gynaecolog Chair P/C	BGY - 401 - PS	
Gynaecolog Chair S/S	BGY -401 - SS	
Gynaecolog Chair With Crank P.C	BGY -40 - CCS	
Gynaecolog Chair With Crank S/S	BGY-201-CC	
Gynaecolog Examination Table P/C + Foot Step	BGT-101	
Gynaecolog Examination Table S/S + Foot Step	BGT - 401 - S	
Gynaecolog Examination Table With Crank P/C + Foot Step	BGT-101-TC	
Gynaecolog Examination Table With Crank S/S	BGT-201-TC	
Gynaecolog Chair Electric	BGY-402-CES	
Gynaecolog Chair BKKBN / Obgyn Bed	BGT-202-BK	
HOSPITAL BED		35
Children Hospital Bed P/C	BHC - 401-S	
Children Hospital Bed S/S	BHC-201	
Children Hospital Bed 1 Crank P/C	BHC-101-1C	
Children Hospital Bed 1 Crank S/S	BHC - 401- 1CR	
Hospital Bed 1 Crank (ABS + Side Rail Alumunium)	BHB-403-1ABSS	
Hospital Bed 1 Crank Deluxe (ABS Panel + Electric)	BHB-101-1E	
Hospital Bed 2 Crank (ABS + Side Rail Alumunium)	BHB-404-2ABSS	
Hospital Bed 2 Crank Elektrik (ABS + Side Rail Alumunium)	BHB-104-2E	
Hospital Bed 2 Crank Hi-Lo	BHB-118-HL	
Hospital Bed 3 Crank (ABS + Side Rail Alumunium)	BHB-405-3ABSS	
Hospital Bed Elektrik	BHB - 408 - ES	
Hospital Verlos Bed P.C + Foot Step	BVB-101	
Hospital Verlos Bed S/S + Foot Step	BVB - 401 - S	
ICU Bed	BIB-101-R	
ICU Bed	BIB-101-E	
Hospital Multifunction Gynaecolog Bed	BBM - 101 - O	

Product Name	Type	Page
INCUBATOR		42
Incubator Digital Servo + Skin Control	BICB-402-SSCS	
Incubator Digital Control	BICB -401- DCS	
Incubator Digital Transport	BICB-104-TR	
INFANT WARMER		44
Infant Warmer + Skin Control Servo	BIW-302	
Infant Warmer Standard	BIW - 401 - S	
INFUSION STAND		45
Infusion Stand 3 Wheel P/C	BIS-101-3W	
Infusion Stand 3 Wheel S/S	BIS-201-3W	
Infusion Stand 4 Wheel P/C	BIS-101-4W	
Infusion Stand 4 Wheel S/S	BIS-201-4W	
Infusion Stand 5 Wheel P/C	BIS - 401 - 5WS	
Infusion Stand 5 Wheel S/S	BIS-201-5W	
Infusion Stand 5 Wheel S/S	BIS-101-5WP	
INSTRUMENT CABINET		47
Instrument Cabinet 1 Pintu P/C	BIC-101-1P	
Instrument Cabinet 1 Pintu S/S	BIC-201-1P	
Instrument Cabinet 2 Pintu P/C	BIC-103-2P	
Instrument Cabinet 2 Pintu S/S	BIC-203-2P	
LEMARI ASAM		48
Lemari Asam (HPL)	BLA-101	
Laminar Air Flow	BLA-101-LAM	
Laminar Air Flow	BLA-101-LAM	
LEMARI NARKOTIKA		
Lemari Narkotika (S) Uk. 20x20x40	BIC-104-N	
Lemari Narkotika (M) Uk. 25x40x60	BIC-105-N	
Lemari Narkotika (L) Uk. 50x60x90	BIC-106-N	

Product Name	Type	Page
MAYO STAND / OVER BED TABLE		49
Mayo Stand S/S	BMS-201	
Over Bed Table P/C	BMS-101	
Over Bed Table Deluxe With Gas Spring P/C	BMS-102-GS	
Mayo Stand S/S	BMS-20-KD	
Over Bed Table P/C	BMS-202	
MEDICINE CABINET		50
Medicine Cabinet 1 Pintu P/C Side Glass	BMC-101-1P	
Medicine Cabinet 1 Pintu S/S	BMC-201-1P	
Medicine Cabinet 2 Pintu P/C	BMC-103-2P	
Medicine Cabinet 2 Pintu S/S	BMC-203-2P	
MORTUARY REFRIGATOR		51
Mortuary Refrigerator 2 pintu	BMO-203-2PR	
MORTUARY CARRIAGE		51
Mortuary Table S/S (Pemandian Jenazah)	BMO-201-T R	
Mortuary Carriage S/S	BMO-202-R	
OPERATING TABLE		52
Operating Table	BOT- 409 - 1ACS	
Operating Table	BOT-210-1AC/DC	
RADIOLOGY TABLE / MEJA X - RAY		53
Radiology Table P/C	BRT - 401 - S	
SCRUB STATION		54
Scrub Station 1 Person (Automatic Sensor)	BSR-401-R	
Scrub Station 2 Person (Automatic Sensor)	BSR-402-2PR	
Scrub Station 3 Person (Automatic Sensor)	BSR-202-3R	
Scrub Station 1 Person (Manual)	BSR-401-1PS	
Scrub Station 2 Person (Manual)	BSR-402-2PS	
Scrub Station 3 Person (Manual)	BSR-202-3PS	
SOUND CHAMBER		56
Sound Proof Chamber Audiometr (Melamic)	BSA-401-M	
Sound Proof Chamber Audiometr (HPL)	BSA-401-T	

Product Name	Type	Page
STAND WASCOM / BOWL STAND		57
Stand Wascom Single P/C	BSW-101-S	
Stand Wascom Single S/S	BSW-201-S	
Stand Wascom Single Deluxe S/S	BSW-201-SD	
Stand Wascom Double P/C	BSW-101-D	
Stand Wascom Double S/S	BSW-201-D	
Stand Wascom Double Deluxe S/S	BSW-201-DD	
DRESSING TROLLEY		58
Dressing Trolley P/C	BDT-101	
Dressing Trolley S/S	BDT-201	
EMERGENCY TROLLEY		58
Emergency Trolley P/C	BETR-101	
Emergency Trolley S/S	BETR-201	
FOOD TROLLEY		59
Dirty Food Trolley	BFT-208-DFT	
USG Trolley Deluxe	BUT-101	
USG Trolley Deluxe	BUT-101-D	
Food Trolley 16 Rack (2 Door) S/S	BFT-202	
Food Trolley 16 Rack (2 Door) S/S Warmer	BFT-202-16S	
Food Trolley 24 Rack (2 Door) S/S	BFT-203	
Food Trolley 30 Rack (2 Door) S/S	BFT-211-30T	
Food Trolley 32 Rack (2 Door) S/S	BFT-204-2F	
INSTRUMENT TROLLEY		62
Instrument Trolley P/C	BIT-102-R	
Instrument Trolley A/S	BIT-202	
Instrument Trolley (Pipa S/S)	BIT-203	
Instrument Trolley 1 Drawer P/C	BIT-102-1DR	
Instrument Trolley 1 Drawer S/S	BIT-202-1DR	
Instrument Trolley 1 Drawer (Laci) P/C	BIT-101-2DR	
Instrument Trolley 2 Drawer S/S	BIT-201-2DR	
Instrument Glass Trolley P/C	BIT-101-R	
Instrument Glass Trolley S/S	BIT-201-R	

Product Name	Type	Page
LAPAROSCOPY TROLLEY		
Laparascopy Trolley	BAT-103-L	
CLOTHES TROLLEY		64
Laundry Cart + Linen P/C	BCT-201	
Laundry Cart + Linen S/S	BCT-202-C	
Laundry Trolley S/S	BLT-210 -T	
Linen Trolley P/C	BLNT-101	
Linen Trolley S/S	BLNT-201	
MEDICINE TROLLEY		66
Medicine Trolley 24 Rak P/C (Acrylic)	BMT-101-24R	
Medicine Trolley 24 Rak S/S (Acrylic)	BMT-201-24R	
Medicine Trolley 24 Rak P/C (Plastik)	BMT-103-24PL	
Medicine Trolley 26 Rak S/S (Stainless Steel)	BMT-201-26R	
Medicine Trolley 30 Rak S/S (Acrylic)	BMT-202-30R	
Medicine Trolley 36 Rak P/C (Plastik)	BMT-103-36PL	
Medicine Trolley 36 Rak P/C (Acrylic)	BMT-101-36R	
OXYGEN TROLLEY		68
Oxygen Trolley 1 m ³ P/C	BOXT-101-CMP	
Oxygen Trolley 1.5 m ³ P/C	BOXT-101-1.5	
Oxygen Trolley 6 m ³ S/S	BOXT-201-6	
Oxygen Trolley 6 m ³ P/C	BOXT-101-6	
U.V MOBILE STERILISATOR		69
U.V Mobile Sterilizer 4 Lamp	BUV-201-4L	
U.V Mobile Sterilizer 4 Lamp Remot	BUV-201-4LR	
U.V Mobile Sterilizer 5 lamp	BUV - 401 - 5LS	
U.V Mobile Sterilizer 5 lamp Remot	BUV-201-5LR	
U.V Safety Cabinet	BUV-202-SC	

Product Name	Type	Page
VISUAL CHART / SNELEN CHART		70
Visual Chart light box / Snelen Chart	BSN-301	
Visual Chart Light Box / Snelen Chart (Remote Control)	BSN-302-R	
X-RAY		70
X-Ray Barrier (Tabir)	BXR-301	
X-RAY VIEWER / LIGHT CASE		71
X-Ray Viewer Single	BXR-101-S	
X-Ray Viewer Single Deluxe	BXR-101-SD	
X-Ray Viewer LED (THIN)	BXR-102-LED	
X-Ray Viewer Double	BXR-101-D	
X-Ray Viewer LED (THIN)	BXR-101-LED	
STANDARD CASSETE		72
Standard Cassete P/C (Knock Down)	BXS-101-KD	

Disinfectant Chamber

BDC - 102 - DC

Swab Chamber

BDC - 103- SC



PRODUCT SPESIFICATION	
Construction	Steel hollow
Material	Acrylic 3 mm, steel hollow 30x30 plywood, & transparent curtain
Dimension	110 (L) x 110 (W) x 240 cm (H)
Tank Capacity	50 L - in (reception), 10 L - out (disposal)
Power Consumption	70 Watt
Note	Automatic sensor, light wirth TL, exhaust fan, nozzle spray, & castor 4"
PRODUCT SPESIFICATION	
Wall Frame	Steel hollow 30x30 mm
Material	Acrylic 5 mm transparent, Shark board 10 mm (PVC food grade)
Dimension	120 (L) x 120 (W) x 220 cm (H)
Power Consumption	70 Watt
Note	Speaker sound 2 way, light with TL, Prefilter , PVC gloves and CCTV

Isolation Chamber

BAS - 702 - IC

Isolation Stretcher

BAS - 502 - IS



PRODUCT SPESIFICATION	
Construction	acrylic 2 mm transparent, Stainless steel hollow 20 x 20, plate t: 1mm,
Mattress	60 mm thickness
Finishing	Polished
Dimension	(L) 210 x (W) 60 x (H) 70 cm,
Power supply	450 watt
Load capacity	Up to 150 Kg
System Filter	Inlet Prefilter and outlet HEPA filter for 0.3um particles
Accessories	Latex gloves, UV lamp, Fan & Battery charger
PRODUCT SPESIFICATION	
Construction	Stainless Steel Pipe ∅ 1"x1.2mm, ∅ 5/8"x1.2mm, Plate L:1mm, t:3mm, Iron As, Aluminium Casting material qualified and selected.
Base	Made of Stainless Steel Plate
Mattress	Mattress foam semi latek covered with vynil, 60 mm thickness
Finishing	Polished
Dimension	(L)190x(W)60x(H)88 cm, (Bendable) (L)190 x(W)60 x(H)88 cm
Power supply	30 Watt-12Volt
Load capacity	Up to 150 Kg
System Filter	Inlet Prefilter and outlet HEPA filter for 0.3um particles
Accessories	Battery and charger

Ambulance & special vehicle

Kami telah berusaha keras dengan yang terbaik dan paling dapat diandalkan untuk menyediakan ambulan dalam industri demi kepuasan pelanggan kami. Pelanggan kami selalu memiliki harapan yang tinggi dan mengharapkan kita untuk bekerja ekstra dalam proses manufaktur dan itulah yang kita lakukan.

Kami juga membangun banyak kendaraan tujuan khusus, di mana pun Anda akan beroperasi, kami akan memastikan bahwa itu dibangun tepat dan sesuai kebutuhan anda.



Spesifikasi Ambulance type Deluxe

- Lampu rotary light bar oval LTF 2000
- Sirene multi suara lengkap dengan mic
- Tempat duduk perawat dilengkapi dengan box peralatan
- Landasan tandu dilengkapi dengan scoop stretcher
- Tabung pemadam kebakaran kapasitas 1 kg
- Lampu sorot belakang model spot light
- 2 buah lampu operasi model geser
- Modifikasi lantai dari plywood dilapisi dengan vinil anti bakteri
- Logo dan tulisan (cutting stiker/digital printing)
- Kaca film 60%-80%
- Lemari peralatan dengan finishing lapis akrilik
- Gantungan infuse model geser terbuat dari bahan stainless steel
- Bumper belakang mobil dilapisi plat stainless steel
- Extra DC
- Central Oxygen BSS system, AM-3001
 - 2 buah tabung oksigen 1m3
 - 2 Tabung reglator standar
 - 1 set selang oksigen tekanan tinggi dengan sistem press
 - 1 set pengatur oksigen dengan kran on/off
 - 1 buah keluaran dinding dengan tulisan "OXIGEN"
- 1 set flowmeter dan humidifier
- 1 buah nasal cannula
- 1 buah masker oksigen

Ambulance stretcher multipurpose BAS - 303 - M, Konstruksi alluminium alloy & Stainless steel complete with mattress,safety belt, infus stand.

- Rail Block With Stoper
- Locking & Unlocking



Spesifikasi Ambulance Type Standard

- * Rotary Flash Light Bar Oval
- * Sirene Multi suara lengkap dengan Mic
- * Tempat duduk perawat dilengkapi dengan Box peralatan
- * Tempat Duduk Perawat / Dokter
- * Landasan Stretcher dilengkapi dengan tempat Scoop Stretcher
- * Tabung pemadam kebakaran kapasitas 0,5 - 1 Kg
- * Lampu sorot belakang model Spot Light
- * 2 (dua) buah lampu Operasi model digeser Halogen
- * Modifikasi lantai dari Plywood dilapisi Vinyl Koin
- * Logo dan tulisan Ambulance Standard dari bahan 3M
- * Kaca Film 60% - 80%
- * Lemari Peralatan dengan Finishing Duco
- * Gantungan Infuse model geser terbuat dari bahan stainless Steel
- * Bumper belakang mobil dilapisi Plat Stainless Steel
- * Extra DC Outlets
- * Oxygen System
 - 1 buah tabung Oxygen 1 m3
 - 1 buah Regulator Standard
 - 1 buah Nassal Cannula
 - 1 buah Masker Oxygen

Ambulance stretcher Standard Folding BAS - 301 - S, Konstruksi alluminium alloy & Stainless steel complete with mattress,safety belt, infus stand.

- Rail Block With Stoper
- Locking & Unlocking



Spesifikasi Ambulance Type Economy

- * Lampu Rotary Light bar Oval
- * Sirene Multi suara lengkap dengan Mic
- * Tempat duduk perawat dilengkapi dengan Box peralatan
- * Tabung pemadam kebakaran kapasitas 0,5 - 1 Kg
- * Lampu sorot belakang model Spot ight
- * Modifikasi lantai dari Plywood dilapisi Vinyl Standard
- * Ambulance teks Standard dari 3M
- * Kaca Film 60% - 80%
- * Gantungan Infuse model geser terbuat dari bahan stainless Steel
- * Kotak P3K + Wastafel Finishing Duco
- * Oxygen System
 - 1 buah tabung Oxygen 1 m3
 - 1 buah Regulator Standard
 - 1 buah Nassal Cannula
 - 1 buah Masker Oxygen
 - Braket Tabung
- Ordinary Stretcher BAS - 201 - S,
Konstruksi Stainless steel complete with mattress,safety belt, infus stand.
- Rail Block With Stoper
- Locking & Unlocking



Spesifikasi Ambulance Jenazah

- * Lampu rotary light oval LTF 2000, AM-6003
- * Sirene multi suara lengkap dengan mic
- * Tempat duduk perawat dilengkapi dengan box peralatan
- * Lampu sorot belakang model spot light
- * Modifikasi lantai dari Plywood dilapisi dengan vinil
- * Logo dan tulisan standar dari bahan stiker 3M
- * Kaca Film 60% - 80%
- * Mortuary Cariatage BMO - 201,
Konstruksi Stainless steel complete with cover and grounding.



Ambulance Stretcher

Stretcher Multipurpose

BAS-403- MS

Specification

Wooden Stretcher
BAS - 307 - K



Folding Stretcher
BAS - 304 - FS



Ordinary Stretcher
BAS - 201 - S



Stretcher Standard
BAS - 302 - S



Specification

Construction : Stainless Steel Hollow,Pipe, Iron As, Plate,Steel,PVC, Aluminium Cor. is a material qualified and selected. (BAS - 302 - S)

Stainless Steel Pipe, Iron As, Plate,PVC, Aluminium Cor. is a material qualified and selected. (BAS - 201 - S & BAS - 304 - FS), Plywood (BAS - 307 - K)

Base Plate : Made of plate Stainless Steel

Mattress : Mattress foam semi latex covered with vinyl, 6cm thickness.

Finishing : Polished

Side Rail : Stainless Steel side rail with handle locks, can be fold and easy it's user

I.V Stand : High-quality Stainless Steel material with 2 Infusion holder

Wheels : RHJTW 4 "using 2 wheels plain and 2 wheel lock

Max Load : 150 Kg

Weight : 35 Kg


Stand size : (L) 190 x (W) 58 x (H) 88 cm (BAS - 302 - S)

(L) 190 x (W) 52 x (H) 36 cm (BAS - 201 - S)

(L) 200 x (W) 60 x (H) 20 cm (BAS - 304 - FS)

(L) 180 x (W) 41 x (H) 4 cm (BAS - 307 - K)

Note: Back Raise,backrest tilt can be adjusted using the regulator is behind the backrest, cast aluminum connection, the legs fitted with PVC to make it easier when the ambulance entrance.



Construction : Hollow pipe, pipe, Plate is a material qualified and selected.

Base Plate : made of Stainless Steel Plate

Mattress : Mattress foam rebonded covered with vinyl, 5cm thickness.

Finishing : Polished

Side Rail : Stainless Steel side rail with handle locks, can be fold and easy its user.

I.V Stand: High-quality Stainless Steel material with 2 Infusion holder

Wheels: 5 "using 2 wheels plain and 2 wheel lock

Max Load: 150 Kg

Weight: 36 Kg

Size : a Sleep : (L) 190 x (W) 58 x (H) 88 cm


Sit Down : (L) 115 x (W) 58 x (H) 140 cm

Note: Position Back Raise, This brancard could be several such positions, folded into a chair with a certain slope, each connection is made of aluminum and iron in Crom, there is an extra wheel to facilitate entry into the ambulance. equipment with bracket oxygen pacric

Stretcher Standard folding

BAS - 401 - SS

Specification



Construction : Pipe, Iron As, PVC, Aluminium Cor, Plate is a material qualified and selected.

Base Plate : Made of plate Stainless Steel

Mattress : Mattress foam rebonded covered with vinyl, 5cm thickness.

Finishing : Polished

Side Rail : Stainless Steel side rail with handle locks, can be fold and easy it's user

I.V Stand : High-quality Stainless Steel material with 2 Infusion holder

Wheels : 5 "using 2 wheels plain and 2 wheel lock

Max Load : 150 Kg

Weight : 34,5 Kg

Stand size : (L) 190 x (W) 60 x (H) 88 cm


Bendable size : (L) 190 x (W) 60 x (H) 36 cm

Note: Back Raise,backrest tilt can be adjusted using the regulator is behind the backrest, cast aluminum connection, the legs fitted with PVC to make it easier when the ambulance entrance. equipment with bracket oxygen pacric


Baby Box



BBB - 101



BBB - 201



Specification

Construction :	Steel Hollow , pipe, Plate is a material qualified and selected for (BBB-101) Stainless Steel Hollow , pipe, Plate is a material qualified and selected for (BBB -201)
Mattress :	Mattress foam semi latex covered with vinyl, 6cm thickness.
Base :	Made of Plywood
Finishing :	Epoxy powder Coating for (BBB-101) Polished for (BBB -201)
Side Rail :	Steel side rail with handle locks, can be fold and easy its user for (BBB-101) Stainless Steel side rail with handle locks, can be fold and easy its user for (BBB-201)
Panel Head & Feet :	Steel for (BBB-201) Stainless Steel for (BBB-201)
Wheels:	Castor 2 "
Max Load:	50 Kg
Size:	(L) 80 x (W) 50 x (H) 120 cm
Note:	On the side rails are 1 piece side access door with a lock so that when the baby is in a safe box.


Baby Cot



BBC - 201



BBC - 401 - S




Specification

Construction :	Steel Hollow , pipe, Plate is a material qualified and selected for (BBC - 401- S) Stainless Steel Hollow , pipe, Plate is a material qualified and selected for (BBC - 201)
Mattress :	Mattress foam semi latex covered with vinyl, 6cm thickness.
Base :	Made of Plywood
Finishing :	Epoxy Powder Coating for (BBC - 401- S) Polished for (BBC - 201)
Side Rail :	Steel side rail with handle locks, can be fold and easy its user for (BBC - 401- S) Stainless Steel side rail with handle locks, can be fold and easy its user for (BBC - 201)
Panel Head & Feet :	Steel for BBC - 401- S Stainless Steel for BBC - 201
Wheels:	Castor 2 "
Max Load:	50 Kg
Size:	(L) 80 x (W) 50 x (H) 120 cm
Note:	On the side rails are 1 piece side access door with a lock so that when the baby is in a safe box.
Accessories :	Infus Stand & Basket for (BBC - 401- S)

Baby Tafel

BTA - 101 - L



BTA - 201



BTA - 401 - LS




BTA - 101



Specification

Construction :	Steel Hollow , Plate, Bar is a material qualified and selected for (BTA - 101 L & BTA - 101) Stainless Steel Hollow , Plate, Bar is a material qualified and selected for (BTA - 401-LS & BTA - 201)
Mattress :	Mattress foam semi latex covered with vinyl, 8cm thickness.
Base Plate :	made of Plywood
Finishing :	Epoxy Powder Coating for (BTA - 101 L & BTA - 101) Plished for (BTA - 401-LS & BTA - 201)
Side Rail :	Steel side rail for (BTA - 101 L & BTA - 101) Stainless teel side rail for (BTA - 401-LS & BTA - 201)
Wheels:	Castor 2 " Castor 3 " (BTA - 401-LS)
Max Load:	50 Kg
Size:	(L) 107 x (W) 58 x (H) 92 cm
Note:	there are 2 pieces for storage drawers, 1cupboard with 2 doors to store clothing and so forth, 1 taray that can be folded when not in use.
Accessories :	1 piece of light for illumination for (BTA - 401-LS & BTA - 101 L)

Baby Measurement Wood



BBT - 302 - K


Specification

Construction :	Plywood is a material qualified and selected.
Base :	Plywood covered with vinyl
Hinge :	Made of Stainless Steel Plate
Finishing :	Varnish
Side Rail :	Made of plywood
size :	Elongated (L) 110 x (W) 33 x (H) 7 cm Folded (L) 57 x (W) 33 x (H) 14 cm
Accessories :	Length gauges baby & lock
Note :	The main ingredient of plywood, there is 1 unit of measure in the baby comes with a limiting, can be folded and fitted lock.
Optional :	-

Baby Measurement Aluminium

BBT - 302 - A


Specification



Construction :	Aluminium, PVC, Acrylic is a material qualified and selected.
Cornerstone :	Made of aluminium
Hinge :	Made of screw
Finishing :	Epoxy powder coating & Screen printing
Delimiter :	Made of Acrylic & PVC
size :	Elongated (L) 105 x (W) 17 x (H) 28 cm Folded (L) 105 x (W) 17 x (H) 5 cm
Accessories :	Rubber cover
Note :	The delimiter can be folded, figures on aluminium using screen printing, at the end there is a rubber delimiter aluminium.
Optional :	-


Bed Screen

BSS - 203 - TR



BSS - 103 - TR

BSS - 202 - DR



BSS - 102 - DR

BSS - 201 - SR



BSS - 101 - SR


Bed Side Cabinet

BSC - 101- A



BSC - 101

BSC - 102 - FT



BSC - 104

BSC - 105 - TO



Specification

Construction :	Steel Pipe,Hollow,Plate is a material qualified and selected.	
Material :	Stainless Steel	(BSS-203-TR / BSS-201-SR / BSS-202-DR)
	Steel	(BSS-103-TR / BSS-101-SR / BSS-102-DR)
Finishing :	Polished	(BSS-203-TR / BSS-201-SR / BSS-202-DR)
	Epoxy Powder Coating	(BSS-103-TR / BSS-101-SR / BSS-102-DR)
Cover :	White cloth	
Type Screen :	Triple	(BSS-203-TR / BSS-103-TR)
	Double	(BSS-102-DR / BSS-202-DR)
	Single	(BSS-101-SR / BSS-201-SR)
Size :	(L) 270 x (W) 50 x (H) 175 cm	(BSS-203-TR / BSS-103-TR)
	(L) 180 x (W) 50 x (H) 175 cm	(BSS-102-DR / BSS-202-DR)
	(L) 90 x (W) 50 x (H) 175 cm	(BSS-101-SR / BSS-201-SR)
Wheels :	2" Wheel, two wheels with locking	
Accessories :	White cloth	
Note :	insulating cloth using cloth adhesive system that can be released and cleaned separately	

Specification

Construction :	Steel Hollow, Iron As, Plate, Plastic ABS is a material qualified.	
Top Cabinet :	Made of plywood and coated in varnish melaminto	(BSC - 101 / BSC - 102 - FT)
	Made of plywood and melaminto	(BSC - 104 / BSC - 105 - TO)
	Made of Plastic ABS	(BSC - 101 - A)
Finishing :	Epoxy Powder Coating & Varnis	
Wheels :	2" Castor, two wheels wearing lockdown.	
Max Load :	50 Kg	
size :	(L) 48 x (W) 46.5 x (H) 81 cm	(BSC - 101 - A)
	(L) 45 x (W) 40 x (H) 90 cm	(BSC - 101)
	(L) 45 x (W) 40 x (H) 95 cm	(BSC - 102 - FT)
	(L) 45 x (W) 40 x (H) 90 cm	(BSC - 104)
	(L) 45 x (W) 40 x (H) 95 cm	(BSC - 105 - TO)
Accessories :	Hanger, Drawer & Food Tray	
Note :	there is 1 piece of drawers, cabinets and 1 lower storage space, bed side there are, hanger towel that can be set	

Blue Light Photo Terapy

Biological Safety Cabinet



Specification

Construction :	Steel Hollow,Plate,Stainless Steel Hollow, Pipe, Iron As,Acrylic is a material qualified
Box Lamp :	Made of Steel Plate
Pole in :	Made of Stainless Steel Hollow
Outer pole :	Made of Steel Hollow
Finishing :	Epoxy Powder Coating
Wheels :	3" Wheel RHJTW, two wheels with locking
Weight :	18 Kg
size :	(L) 70 x (W) 60 x (H) 120 - 180 cm
Accessories :	Timer, Switch on / off, Ballast Electric,TL Blue Philip.
Lamp :	20 Watt 1 Lamp
Power Consumption :	80 Watt, 4 x 20 watt (BPT - 301 - 4L) 100 Watt, 5 x 20 Watt (BPT - 301 - 5L)
Power Supply :	AC 220 Volt, 50/60 Hz
Note :	On the part of the pole there is a retaining spring. in any moving parts fitted lock.

BSF - 101 - C

BSF - 101



Specification


Construction :	Steel Hollow ,Stainless Steel Plate, Steel Plate,Acrylic, is a material qualified and selected.
Base :	Stainless Steel Plate 304 t:1mm
Finishing :	Epoxy Powder Coating & Polished (BSF - 101 - C) Epoxy Powder Coating (BSF - 101)
Door :	slide up and comes with transparent acrylic
Handle :	Aluminium
Max Load :	50 Kg
size :	(L) 136 x (W) 80 x (H) 155 cm (BSF - 101 - C) (L) 62 x (W) 52 x (H) 140 cm (BSF - 101)
Accessories :	Door locks and handles, blower, U.V lamp, Lighting
Power Consumption :	450 Watt (BSF - 101 - C) 24 Watt (BSF - 101)
Input Voltage :	220 Volt, 50/60 Hz
Hepa Filter	Double Hepa Filter 0.3 um for (BSF - 101 - C)
Note :	The exterior is made of steel plate with epoxy powder coating finishing, and for the interior are made of stainless teel plate with polished finishing.

Brangcart



BBP - 403 - CS

BBP - 103 - C



Specification

Construction :	Hollow, pipe,Iron Bar, Plate is a material qualified and selected for
Material :	Stainless Steel (BBP - 403 - CS) Steel (BBP - 103 - C)
Matrress :	Matrress super foam covered with vinyl, 8cm thickness.
Base :	made of Plywood
Finishing :	Polished for (BBP - 403 - CS) Epoxy Powder Coating (BBP - 103 - C)
Side Rail :	Stainless Steel side rail with locks for BBP - 403 - CS Steel side rail with locks for BBP - 103 - C
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels:	RHJTW 4 "using 2 wheels plain and 2 wheel lock (BBP - 103 - C) RHJTW 5 "using 2 wheels plain and 2 wheel lock (BBP - 403 - CS)
Position :	Back Raise with manual adjuster (BBP - 103 - C) Back Raise with Gaspring adjuster (BBP - 403 - CS)
Max Load:	150 Kg
Size:	(L) 200 x (W) 77 x (H) 80 cm
Note:	section side rail can be folded or lowered with the manual system and the tilt adjustable backrest section with backraise adjuster.

Brangcart



BBP - 102 - S




BBP - 101



BBP - 401 - S



BBP - 202 - S



Specification

Construction :	Hollow, pipe,Iron Bar, Plate is a material qualified and selected for
Material :	Stainless Steel (BBP - 401 - S & BBP - 202 - S) Steel (BBP - 101 & BBP - 102 - S)
Matrress :	Matrress Super foam covered with vinyl, 8cm thickness.
Base :	made of Plywood
Finishing :	Polished for (BBP - 401 - S & BBP - 202 - S) Epoxy Powder Coating (BBP - 101 & BBP - 102 - S)
Side Rail :	Stainless Steel side rail with locks for (BBP - 401 - S & BBP - 202 - S) Steel side rail with locks for (BBP - 101 & BBP - 102 - S)
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels:	RHJTW 4 "using 2 wheels plain and 2 wheel lock
Max Load :	150 Kg
Size :	(L) 200 x (W) 77 x (H) 80 cm
Note :	section side rail can be folded or lowered with the manual system and the tilt adjustable backrest section with backraise adjuster.
Other Excess :	brangcart This can also be enabled for the stretcher, brangcart body parts and a separate pedestal for (BBP - 102 - S & BBP - 202 - S)
Accessories :	Basket for (BBP - 401 - S)

Brangcart

Phlebotomi Chair

Brangcart Complete Hig-low Deluxe
BBP - 111 - DA



Transper Strecher Trolley
BBP - 110 - TST



Brangcart Hig-low Deluxe
BBP - 104 - HLR



BPC - 101



BPC - 201



Phlebotomi Chair Electric
BPC - 201 - E



Specification



Construction :	Steel Hollow pipe, pipe, Iron bar, Stainless Steel Pipe, Plate is a material qualified and selected.
Base Plate :	made of sheet plate (SPCC) the perforated with embossing process (BBP-104-D, BBP-110-TST) made of ABS Molding (BBP-111-DA)
Finishing :	Epoxy Powder Coating & Polished
Side Rail :	Big ABS panel side rail single with locks, can be fold and easy its user. (BBP-111-DA) Small ABS panel side rail single with locks, can be fold and easy its user. (BBP-104-D) Aluminium side rail single with locks, can be fold and easy its user. (BBP-110-TST)
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels :	4 "using 2 wheels plain and 2 wheel lock
Movers :	Gastpring & Manual Crank Adjuster
Max Load :	150 Kg
Size :	(L) 190 x (W) 66 x (H) 51 - 85 cm (BBP-111-DA) (L) 200 x (W) 65 x (H) 55 - 87 cm (BBP-104-D) (L) 200 x (W) 65 x (H) 55 - 87 cm (BBP-110-TST)
Accessories :	Mattress, Infus Stand, Bracket Oxygen (BBP-104-D) Mattress, Infus Stand, Bracket Oxygen, Basin made of PVC, Bumper protectors (BBP-110-TST & BBP-111-DA)
Note :	Back Raise, High & Low, This brangcart can be height adjustable with a crank, the tilt adjustable backrest section with backraise Gastpring. There are 1 piece bracket oxygen, at the bottom there is a body made of PVC as well as the container.

Specification

Construction :	Steel Hollow, Pipe, plate, Stainless Steel Pipe, PVC, Iron As, (BPC-201-E) Hollow, Pipe, plate, Iron As, is a material qualified and selected. (BPC-101 / BPC-201)
Material :	Stainless steel (BPC-201) Steel (BPC-101 / BPC-201-E)
Mattress :	Mattress foam semi latex covered with vinyl.
Retaining Arms :	Frame made of Stainless Steel
Finishing :	Polished (BPC-201) Epoxy Powder Coating (BPC-101 / BPC-201-E)
Footwear :	Made from a rubber
Max Load :	150 Kg
Size :	(L) 75 x (W) 75 x (H) 115 cm (BPC-201-E) (L) 95 x (W) 60 x (H) 120 cm (BPC-101 / BPC-201)
Movers :	Actuator 4000 N / 24V, 3.8 A (BPC-201-E)
Movement Control :	Hand Controller
Note :	the slope of the rear backrest and footrest can be adjusted front, the sides are supporting hand side barrier as well. (BPC-201-E) the slope of the rear backrest and footrest can be adjusted front, the sides are supporting hand side barrier as well. (BPC-101 / BPC-201)

Wheel Chair

BKR - 203 - PP



BKR - 402 - S



Specification

Construction :	Pipe, Plate is a material qualified and selected.
Base :	Syntek Leater
cushion :	Synthetic leather is supple and strong
Finishing :	Cromed (BKR - 402-S) Epoxy Powder Coating (BKR - 203)
Side Rail :	Steel Pipe and plastic
Mover :	Wheel
Rim :	Fingers (BKR - 402-S) Plastic Mold (BKR - 203)
Wheels :	Rear wheels 21" and front wheels 8"
Max Load :	100 Kg
size :	(L) 85 x (W) 63 x (H) 90 cm
Accessories :	Rubber Wheel, custody feet.
custody feet :	Can be folded or removed
Note :	On the rear wheel brakes on both wheels fitted, this wheelchair can be folded when not in use.
Accessories :	Infus Stand + identity boards (BKR - 402-S)

Examination Table

BET - 102 - A




BET - 401 - PS



BET - 401 - SS



BET - 201 - E



Specification


Construction :	Stainless Steel Hollow pipe, pipe, iron bar ,Plate is a material qualified and selected.
Material :	Stainless Steel (BET - 401-SS & BET - 201- E) Steel (BET - 401-PS & BET - 102 - A)
Finishing :	Polished (BET - 401-SS & BET - 201- E) Epoxy Powder Coating (BET - 401-PS & BET - 102 - A)
Mattress :	Mattress Super foam covered with vinyl, 8cm thickness.
Base Plate :	made of Plywood
Footwear :	Made from a rubber mold
Movers :	Actuator Motor (BET - 201- E) Manual Adjuster (BET - 401-SS,BET - 401-PS,BET - 102 - A)
Max Load :	150 Kg
Size :	(L) 190 x (W) 65 x (H) 80 cm (BET - 401-SS,BET - 401-PS,BET - 102 - A) (L) 200 x (W) 65 x (H) 55 - 87 cm (BET - 201- E)
Note :	Position Back Raise,High & Low, The examination table can be height adjustable and backrest sections can also be set by using the tilt actuator propulsion engines with a power 4000 - 6000 N. (BET - 201- E)
Position :	Back Raise (BET - 401-SS,BET - 401-PS,BET - 102 - A)
Accessories :	There are 3 small drawers and 1 large (BET - 102 - A)
	Include Foot Step (BET - 401-SS,BET - 401-PS,BET - 102 - A, BET - 201- E)
	Small Drawer (BET - 401-SS,BET - 401-PS)

Examine Lamp LED


BEL - 104 - LD



BEL - 102 - LED




BEL - 102 - H



Specification

Construction :	Steel Pipe, Plate is a material qualified and selected. (BEL-104-LD)
	Stainless Steel Pipe, Plate is a material qualified and selected. (BEL-102-LED & BEL-102-H)
Foot :	Made of Steel pipe (BEL-104-LD)
	Made of Plastic (BEL-102-LED & BEL-102-H)
Finishing :	Epoxy Powder Coating (BEL-104-LD)
	Polished (BEL-102-LED & BEL-102-H)
Pole :	Steel pipe (BEL-104-LD)
	Stainless steel pipe (BEL-102-LED & BEL-102-H)
Wheels :	2", 2 wheels wearing locks
size :	(L) 56x (W) 51 x (H) 214 cm (BEL-104-LD)
	(L) 55 x (W) 55 x (H) 209 cm (BEL-102-LED & BEL-102-H)
Power Consumption :	5 Watt (BEL-102-LED & BEL-102-H)
	25 watt (BEL-104-LD)
Input voltage :	AC 220 - 230 Volt, 50/60 Hz
out put :	AC 12 Volt, 5A
Accessories :	LED lamp (BEL-102-LED, BEL-104-LD)
	Halogen Lamp (BEL-102-H)
Note:	pole and lamp system uses a knock down, the neck of the lamp is made of a material that is strong and flexible for ease of movement.


Foot Step P/C



Specification

Construction :	Hollow, Plywood is a material qualified and selected. (BFS-101)
	pipe, Plywood is a material qualified and selected. (BFS-201)
Material :	Steel (BFS-101)
	Stainless Steel (BFS-201)
Base Plate :	Made of Plywood covered with vinyl
Finishing :	Epoxy Powder Coating (BFS-101)
	Polished (BFS-201)
Max Load :	100 Kg
size :	(L) 50 x (W) 40 x (H) 35 cm
Accessories :	Pedestal footing and plastic shoe ladder
Note :	the bottom of the foot there is a protective plastic shoes, footing pads are made of plywood coated with vinyl


Gynaecolog Chair Electric




BGY - 402 - CES

Specification

Construction :	Stainless Steel Pipe, Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Plywood
Matress :	Matress foam semi latex covered with vinyl, 6cm thickness.
Finishing :	Polished
Side Rail :	Stainless Steel side rail
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Actuator Motor
Max Load :	150 Kg
Weight :	76 Kg
size :	(L) 190 x (W) 67 x (H) 67 - 87 cm
Accessories :	Foot Rest, Leg Support, Foot Step, Side Guard Stainless Steel, Dirt Basin.
Position :	Back Raise, Seat, High & Low,
Note:	Parts of the junction legs/foot removable backrest, Actuator with power 4000N-6000N, Ac 110/220 Volt 50/60Hz & 300 VA. Prisoners thigh and handrails , accessories can be shifted

Gynaecolog Chair
BGY - 401 - CCS

BGY - 401 - PS




BGY - 101 - CC

BGY - 401 - SS

BGT - 202 - BK

Specification

Construction	Hollow , Pipe,Iron bar,Plate, plate elbow is a material qualified and selected.
Material	Steel (BGY-401-PS, BGY - 101 - CC) Stainless Steel (BGY-401-SS, BGY - 401 - CCS, BGT - 202 - BK)
Mattress	Mattress foam semi latex covered with vinyl, 8cm thickness.
Base	made of Plywood
Finishing	Epoxy Powder Coating (BGY-401-PS, BGY - 101 - CC) Polished (BGY-401-SS, BGY - 401 - CCS, BGT - 202 - BK)
Footwear	Made from a rubber mold
Dimensions	(L) 170 x (W) 65 x (H) 90 cm (BGY-401-PS, BGY-401-SS, BGY - 401 - CCS,BGY - 101 - CC) (L) 165 x (W) 65 x (H) 80 cm (BGT - 202 - BK)
Max. Load	150 Kg
Side Guard	pipe side rail
Mover	By manual Adjuster & Gaspring (BGY-401-PS, BGY-401-SS) By manual Adjuster ,Crank & Gaspring (BGY - 401 - CCS,BGY - 101 - CC)
Position	Backrest, Knee raise (BGY-401-PS, BGY-401-SS, BGT - 202 - BK) Backrest, Knee raise & Cushion raise (BGY - 101 - CC, BGY - 401 - CCS)
Note	the backrest and footrest can be adjusted by using the tilt adjuster ,the adjuster is made of SPHC plate with a thickness of 5 mm and miling process, To ease movement / tilt, cushion using crank.
Accessories	Leg Support, Holder the basin, Stand Infus,


Gynaecolog Examine
BGT - 401 - S

BGT - 101

BGT - 101 - TC

BGT - 201 - TC

Specification


Construction	Hollow , Pipe , pipe ,Plate , Iron bar material qualified and selected.
Material	Stainless Steel (BGT - 401 - S & BGT - 201 - TC) Steel (BGT - 101 & BGT - 101 - TC)
Base	Made of Plywood
Mattress	Mattress foam semi latek covered with vinyl, 80 mm thickness
Finishing	Epoxy Powder Coating (BGT - 101 & BGT - 101 - TC) Polised (BGT - 401 - S & BGT - 201 - TC)
Dimensions	L. 1800 x W. 650 x H. 880 mm
Max. Load	150 Kg
Infus Stand	High quality Stainless Steel material with 2 infused holder
Mover	By manual Adjuster (BGT - 101) By manual Crank & Adjuster (BGT - 201 - TC & BGT - 101 - TC) By manual Adjuster & Gaspring (BGT - 401 - S)
Position	Backrest, Head raise (BGT - 101 & BGT - 401 - S) Backrest, Head raise & Cushion raise (BGT - 201 - TC & BGT - 101 - TC)
Note	The backrest and Head raise can be adjusted by using the tilt adjuster ,the adjuster is made of SPHC plate with a thickness of 5 mm , To ease movement / tilt, cushion using crank.
Accessories	Leg Support, Holder the basin, Stand Infus, Foot Step

Hospital Bed Children


BHC - 101 - 1C




BHC - 401 - 1CR



BHC - 401 - S



BHC - 201




Specification

Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Material:	Steel (BHC - 401 - S,BHC-101-1C) Stainless Steel (BHC-201,BHC - 401 - 1CR)
Base Plate :	Made of Steel Sheet
Mattress :	Mattress foam semi latex covered with vinyl, 11cm thickness.
Finishing :	Epoxy Powder Coating (BHC - 401 - S,BHC-101-1C) Polished (BHC-201,BHC - 401 - 1CS)
Side Rail :	Pipe
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Manual Crank (BHC-101-1C,BHC - 401 - 1CS)
Wheels :	4" RHJTW, two-wheel and 2 wheel locks wear plain 4" two-wheel and 2 wheel locks wear plain, anti hair wheel (BHC - 401 - S,BHC - 401 - 1CS)
Max Load :	100 Kg
size :	(L) 180 x (W) 82 x (H) 135 cm
Accessories :	Side guard, Stand Infus
Foot & Head Panel :	Made of acrylic with character drawing (BHC - 401 - S,BHC - 401 - 1CS)
Note:	the base bed using knockdown system, to raise and lower side wear guard manual system. backrest is moved by using a crank

Hospital Bed 1 Crank ABS

BHB - 403 - 1ABSS

Specification



Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Steel Sheet
Mattress :	Mattress foam semi latex covered with vinyl, 12cm thickness.
Finishing :	Epoxy Powder Coating
Side Guard :	Aluminium
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Manual Crank made of Plastic
Wheels :	5" two-wheel and 2 wheel locks wear plain
Max Load :	150 Kg
Weight :	60 Kg
size :	(L) 220 x (W) 95 x (H) 55 cm
Accessories :	Aluminium Side guard, Stand Infus.
Position :	Backrest
Note:	the side guard is equipped with a lock, head and foot panels use knockdown system.

Hospital Bed 1 Crank Electric

BHB - 101 - 1E

Specification



Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Steel Sheet
Mattress :	Mattress foam semi latex covered with vinyl, 12cm thickness.
Finishing :	Epoxy Powder Coating
Side Guard :	Aluminium
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	Electrical Actuator Motor
Wheels :	4" RHJTW, two-wheel and 2 wheel locks wear plain
Max Load :	150 Kg
Weight :	62 Kg
size :	(L) 220 x (W) 95 x (H) 55 cm
Accessories :	Aluminium Side guard, Stand Infus.
Position :	Backrest
Note:	the side guard is equipped with a lock, head and foot panels use knockdown system.

Hospital Bed 2 Crank ABS

BHB - 404 - 2ABSS

Specification



Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Steel Sheet
Mattress :	Mattress foam semi latex covered with vinyl, 12cm thickness.
Finishing :	Epoxy Powder Coating
Side Guard :	Aluminium
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Manual Crank Made of Plastic
Wheels :	5" two-wheel and 2 wheel locks wear plain anti hair wheel
Max Load :	150 Kg
Weight :	74 Kg
size :	(L) 220 x (W) 95 x (H) 60 cm
Accessories :	Aluminium Side guard, Stand Infus.
Position :	Backrest, Knee Raise
Note:	the side guard is equipped with a lock, head and foot panels use knockdown system.




Hospital Bed 2 Crank Hi Low

BHB - 118 - HL

Specification



Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Steel Sheet
Mattress :	Mattress foam semi latex covered with vinyl, 12cm thickness.
Finishing :	Epoxy Powder Coating
Side Guard :	Aluminium
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Manual Crank
Wheels :	4" RHJTW, two-wheel and 2 wheel locks wear plain
Max Load :	150 Kg
Weight :	76 Kg
size :	(L) 220 x (W) 95 x (H) 46 - 68 cm
Accessories :	Aluminium Side guard, Stand Infus.
Position :	Backrest, High Low
Note:	the side guard is equipped with a lock, head and foot panels use knockdown system.



Hospital Bed 2 Crank Electric

BHB - 104 - 2E

Specification



Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Steel Sheet
Mattress :	Mattress foam semi latex covered with vinyl, 12cm thickness.
Finishing :	Epoxy Powder Coating
Side Guard :	Aluminium
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	Electrical Actuator Motor
Wheels :	4" RHJTW, two-wheel and 2 wheel locks wear plain
Max Load :	150 Kg
Weight :	62 Kg
size :	(L) 220 x (W) 95 x (H) 60 cm
Accessories :	Aluminium Side guard, Stand Infus.
Position :	Backrest, Knee Raise
Note:	the side guard is equipped with a lock, head and foot panels use knockdown system.




Hospital Bed 3 Crank ABS

BHB - 405 - 3ABSS

Specification



Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Steel Sheet
Mattress :	Mattress foam semi latex covered with vinyl, 12cm thickness.
Finishing :	Epoxy Powder Coating
Side Guard :	Aluminium
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Manual Crank Made of Plastic
Wheels :	5" two-wheel and 2 wheel locks wear plain anti hair wheel
Max Load :	150 Kg
Weight :	76.5 Kg
size :	(L) 220 x (W) 95 x (H) 46 - 68 cm
Accessories :	Aluminium Side guard, Stand Infus.
Position :	Backrest, Knee Raise & High Low
Note:	the side guard is equipped with a lock, head and foot panels use knockdown system.




Hospital Bed Electric

BHB - 408 - ES

Specification

Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Steel Sheet
Mattress :	Mattress foam semi latex covered with vinyl, 12cm thickness.
Finishing :	Epoxy Powder Coating
Side Guard :	Aluminium
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Electrical Actuator Motor
Wheels :	5" two-wheel and 2 wheel locks wear plain anti hair wheel
Max Load :	200 Kg
Weight :	108 Kg
size :	(L) 220 x (W) 95 x (H) 55 - 87 cm
Accessories :	Aluminium Side guard, Stand Infus.Bumper
Position :	Backrest, Knee Raise & High Low
Note :	the side guard is equipped with a lock, head and foot panels use knockdown system.
Optional :	Electrical use remote control




Hospital Verlos Bed

BVB - 401 - S

Specification

Construction :	Stainless Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Plywood
Mattress :	Mattress foam semi latex covered with vinyl, 6cm thickness.
Finishing :	Polished
Side Rail :	Stainless Steel side rail
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels :	4" two-wheel and 2 wheel locks wear plain
Max Load :	150 Kg
Weight :	36 Kg
size :	(L) 208 x (W) 83 x (H) 120 cm
Accessories :	Leg Support, Holder the basin, Stand Infus, Foot Step, Bracket Oxygen
Note :	the backrest can be adjusted by using the tilt adjuster, the adjuster is made of SPHC plate with a thickness of 5 mm and milling process.
Optional :	The bed can be divided into two parts.
Optional :	all parts of the system uses a knock down.



Hospital Verlos Bed

BVB - 101

Specification

Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of Plywood
Mattress :	Mattress foam semi latex covered with vinyl, 6cm thickness.
Finishing :	Epoxy Powder Coating
Side Rail :	Steel side rail
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels :	2" wheel
Max Load :	150 Kg
Weight :	36.5 Kg
size :	(L) 208 x (W) 83 x (H) 120 cm
Accessories :	Leg Support, Holder the basin, Stand Infus, Foot Step.
Note :	the backrest can be adjusted by using the tilt adjuster, the adjuster is made of SPHC plate with a thickness of 5 mm and milling process.
Optional :	The bed can be divided into two parts.
Optional :	all parts of the system uses a knock down.




ICU Bed

BIB - 101 - R

Specification

Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of ABS
Mattress :	Mattress foam semi latex covered with vinyl, 12cm thickness.
Finishing :	Epoxy Powder Coating
Side Guard :	ABS Side Guard
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Electrical Actuator Motor
Wheels :	4" wheel and Central locks
Max Load :	200 Kg
Weight :	108 Kg
size :	(L) 220 x (W) 95 x (H) 55 - 87 cm
Accessories :	ABS Side guard, Stand Infus.
Position :	Backrest, Knee Raise & High Low
Note :	the side guard is equipped with a lock, head and foot panels use knockdown system.
Optional :	Electrical use remote control




ICU Bed Deluxe

BIB - 101 - E

Specification



Construction :	Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.
Base Plate :	Made of ABS
Mattress :	Mattress foam semi latex covered with vinyl
Finishing :	Epoxy Powder Coating
Side Guard :	ABS Side Guard
Head and Foot Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Electrical Actuator Motor
Wheels :	4" Suppo, two-wheel and 2 wheel locks wear plain
Max Load :	200 Kg
Weight :	110 Kg
size -	(L) 220 x (W) 95 x (H) 55 - 87 cm
Accessories :	ABS Side guard, Stand Infus.
Position :	Backrest, Knee Raise,High Low, Trendeleunbrug, Anti Trendeleunbrug
Note:	the side guard is equipped with a lock, head and foot panels use knockdown system. Electrical use remote control




Multifunction Gynaecolog Bed

BBM - 101 - O

Specification



Construction :	Steel Pipe,Hollow, Iron bar, Plate is a material qualified and selected.
Base :	Made of plate
Mattress :	Mattress foam semi latex covered with vinyl
Finishing :	Epoxy Powder Coating
Side Guard :	ABS Side Guard
Head Panel :	ABS Panel
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Mover :	By Electrical Actuator Motor
Wheels :	4" Suppo, two-wheel and 2 wheel locks wear plain
Max Load :	200 Kg
Size :	(L) 200 x (W) 95 x (H) 68 - 92 cm
Accessories :	ABS Side guard, Stand Infus, Basin, Retaining thigh,Trolley storage components.
Position :	Backrest, Knee Raise,High Low, Trendeleunbrug, Anti Trendeleunbrug
Note:	Foot section bearing components can be removed when not in use



Incubator Servo Control

BICB - 402 - SSCS



Specification

Construction :	Steel Hollow, Plate, Stainless Steel Pipe,Hollow, Acrylic is a material qualified and selected.
Base Plate :	Made of Acrylic & Plate
Mattress :	Mattress foam semi latex covered with vlanel, 8 cm thickness.
Finishing :	Epoxy Powder Coating & polished
Side Rail :	Acrylic 3mm
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels :	2" Wheel
Max Load :	8 Kg
Weight :	80 Kg
size :	(L) 110 x (W) 64 x (H) 143 cm
Accessories :	Thermometer, Higrometer, Stand Infus,Mattress, Bracket Oxigen.
Power Consumption :	350 Watt, AC 220 Volt, 50/60 Hz
Temperature Control:	25°c - 38°c Elektronik Profesional Control High Air Temperature 38°c Audible with buzer
Note :	equipped with an alarm system and servo control skin (skin sensor), equipped with iris


Incubator Digital Control

BICB - 401 - DCS

Specification



Construction :	Steel Hollow, Plate, Stainless Steel Pipe,Hollow, Acrylic is a material qualified and selected.
Base Plate :	Made of Acrylic & Plate
Mattress :	Mattress foam semi latex covered with vlanel, 8 cm thickness.
Finishing :	Epoxy Powder Coating & polished
Side Rail :	Acrylic 3mm
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels :	4" Wheel Suppo, two wheels with locking
Max Load :	8 Kg
Weight :	70 Kg
size :	(L) 105 x (W) 51 x (H) 140 cm
Accessories :	Thermometer, Higmeter, Stand Infus,Mattress, Bracket Oxigen.
Power Consumtion :	350 Watt
Power Supply :	AC 220 Volt, 50/60 Hz
Temperature Control:	25 ^o c - 38 ^o c Elektronik Profesional Control High Air Temperature 38 ^o c Audible with buzzer
Note :	equipped with an alarm system




Infant Warmer Digital Servo

BIW - 302

Specification



Construction :	Steel Hollow,Plate,Stainless Steel Plate, Pipe ,Iron As,Acrylic, PVC, is a material qualified Made of Steel Plate
Base Plate :	Mattress foam semi latex covered
Mattress :	with vlanel, 6 cm thickness. Epoxy Powder Coating
Finishing :	Stainless Steel & Acrylic 3mm
Side Rail :	High-quality Stainless Steel material with
I.V Stand :	2 Infusion holder
Wheels :	4" Wheel RHJTW, two wheels with locking
Max Load :	10 Kg
Weight :	50 Kg
size :	(L) 70 x (W) 60 x (H) 180 cm
Accessories :	Thermometer, Higmeter, Stand Infus, Mattress, Cabinet.
Power Consumtion :	600 Watt
Power Supply :	AC 220 Volt, 50/60 Hz
Temperature Control:	25 ^o c - 38 ^o c Elektronik Profesional Control High Air Temperature 38 ^o c Audible with buzzer
Note :	equipped with an alarm system and Servo skin control (skin sensor)




Incubator Digital Transport

BICB - 104 - TR

Specification



Construction :	Steel Hollow, Pipe,Plate,Acrylic is a material qualified and selected.
Base Plate :	Made of Acrylic & Plate
Mattress :	Mattress foam semi latex covered with vlanel, 8 cm thickness.
Finishing :	Epoxy Powder Coating
Side Rail :	Acrylic 3mm
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels :	4" Wheel Suppo, two wheels with locking
Max Load :	8 Kg
Weight :	57.5 Kg
size :	(L) 89 x (W) 45 x (H) 113 - 144 cm
Accessories :	Thermometer, Higmeter, Stand Infus, Mattress, Accu 12 Volt.
Power Consumtion :	350 Watt
Power Supply :	AC 220 Volt, 50/60 Hz
Temperature Control:	25 ^o c - 38 ^o c Elektronik Profesional Control High Air Temperature 38 ^o c Audible with buzzer
Note :	equipped with an alarm system.




Infant Warmer Digital Standard

BIW - 401 - S


Specification



Construction :	Steel Hollow,Plate,Stainless Steel Plate, Pipe ,Iron As,Acrylic, PVC, is a material qualified and selected.
Base Plate :	Made of Steel Plate
Mattress :	Mattress foam semi latex covered with vlanel, 6 cm thickness.
Finishing :	Epoxy Powder Coating
Side Rail :	Stainless Steel & Acrylic 3mm
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder
Wheels :	4" Wheel RHJTW, two wheels with locking
Max Load :	10 Kg
Weight :	53 Kg
size :	(L) 70 x (W) 60 x (H) 180 cm
Accessories :	Thermometer, Higmeter, Stand Infus, Mattress, Cabinet. action lights
Power Consumtion :	600 Watt
Power Supply :	AC 220 Volt, 50/60 Hz
Temperature Control:	25 ^o c - 38 ^o c Elektronik Profesional Control High Air Temperature 38 ^o c Audible with buzzer
Note :	equipped with an alarm system.




Infusion Stand



Specification

Construction :	Steel Pipe,Hollow. Plate,Stainless Steel Pipe, Iron As is a material qualified and selected.
Material :	Steel
Foot :	Made of Steel Hollow (BIS-101-3W,BIS-101-4W,BIS-401-5WS) Made of Plastic molded (BIS-101-5WP)
Finishing :	Epoxy Powder Coating & Polished
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder (BIS-101-3W,BIS-101-4W) High-quality Stainless Steel material with 4 Infusion holder (BIS-401-5WS)
Mover :	By Manual
Wheels :	2" Castor, two wheels wearing lockdown.
Max Load :	2 Kg
size :	(L) 48 x (W) 36 x (H) 112 - 200 cm (BIS-101-3W) (L) 50 x (W) 50 x (H) 112 - 200 cm (BIS-101-4W) (L) 50 x (W) 50 x (H) 112 - 194 cm (BIS-401-5WS)
Accessories :	Locking pole
Note:	the pole in the infusion can be height adjustable up to 190 and use the locking knobstar

Infusion Stand



Specification

Construction :	Stainless Steel Pipe,Plate, Iron As is a material qualified and selected.
Material :	Stainless Steel
Foot :	Made of Stainless Steel Pipe Made of Plastic molded
Finishing :	Polished
I.V Stand :	High-quality Stainless Steel material with 4 Infusion holder
Mover :	By Manual
Wheels :	2" Castor, two wheels wearing lockdown.
Max Load :	2 Kg
size :	(L) 48 x (W) 36 x (H) 112 - 200 cm (BIS-201-3W) (L) 50 x (W) 50 x (H) 112 - 200 cm (BIS-201-4W) (L) 50 x (W) 50 x (H) 112 - 194 cm (BIS-201-5W & BIS-201-5WP)
Accessories :	Locking pole
Note :	the pole in the infusion can be height adjustable up to 190 and use the locking knobstar

Instrumen Cabinet



Instrument Cabinet 1 Pintu
BIC - 201 - 1P



Instrument Cabinet 2 Pintu
BIC - 103 - 2P



Instrument Cabinet 2 Pintu
BIC - 203 - 2P



Instrument Cabinet 1 Pintu
BIC - 101 - 1P

BIC - 102 - 1P BIC - 101 - 1P


Specification

Construction :	Stainless Steel Hollow, Plate is a material qualified and selected.	
Material :	Steel	(BIC-101-1P, BIC-102-1P & BIC-103-2P)
	Stainless Steel	(BIC-201-1P & BIC-203-2P)
Shelves :	4 Shelves Made of Glass	
Braket :	Made of Stainless Steel, retention bracket height can be adjusted according to need	
Type door :	Double Door	(BIC-203-2P & BIC-103-2P)
	Single Door	(BIC-101-1P, BIC-102-1P & BIC-201-1P)
Finishing :	Epoxy Powder Coating	(BIC-101-1P, BIC-102-1P & BIC-103-2P)
	Polished	(BIC-201-1P & BIC-203-2P)
Side Guard :	Made of Glass	(BIC-101-1P)
	Made of plate	(BIC-201-1P, BIC-203-2P, BIC-102-1P & BIC-103-2P)
Max Load :	20 Kg	
size :	(L) 70 x (W) 45 x (H) 170 cm	(BIC-101-1P , BIC-201-1P, BIC-102-1P)
	(L) 90 x (W) 45 x (H) 170 cm	(BIC-203-2P & BIC-103-2P)
Accessories :	Lock door, Handle, Bracket Stainless Steel.	
Note:	the side of the cabinets are made of Glass, shelves made of glass 5mm and adjustable, equipped with lock and handle.	

Lemari Asam (Fume Hood)



BLA - 101 BLA - 101 - LAM




Specification

Construction :	Plywood, HPL, Acrylic, Melamin is a material qualified and selected.
Base :	Made of Plywood Covered HPL
Finishing :	HPL
Door :	slide up and comes with transparent acrylic
Handle :	Aluminium
Max Load :	25 Kg
size :	(L) 120 x (W) 80 x (H) 240 cm
Accessories :	Wastafel, 1 bottom cabinet, door locks and handles, blower, U.V lamp
Power Consumption :	450 Watt
Input Voltage :	220 Volt, 3A, 50/60 Hz
Note :	the exterior is made of plywood HPL coated, and the inside is coated with melamine, there is a sink, lighting uses UV light.
System :	Knock down (BLA - 102 -KD)

Lemari Narkotika


BIC - 104 - N

Specification



Construction :	Steel Plate is a material qualified and selected.
Base :	Made of Steel Plate
Finishing :	Epoxy Powder Coating & Polished
Foot :	Rubber
Max Load :	25 Kg
Size :	BIC-104-N "S" - (L) 20 x (W) 20 x (H) 40 cm BIC-105-N "M" - (L) 25 x (W) 40 x (H) 60 cm BIC-106-N "L" - (L) 50 x (W) 60 x (H) 90 cm
Accessories :	Handle, Lock door
Door :	There are two security doors fitted with locks on every door.
Selves :	There are 3 Selves, 4 storage space.
Note :	has 2 doors and locks, 3-shelf parts made of rubber feet.
Optimal :	Can adjust to customer needs all materials are made of stainless steel


Over Bed Table / Mayo Stand



BMS - 201- KD



BMS - 101




BMS - 201

Over Bed Table Deluxe
BMS - 102 - GS



BMS - 101 - ABS



Specification

Construction :	Pipe, Hollow, iron bar, plate, Plywood & melaminto is a material qualified and selected.	
Material :	Stainless Steel	(BSM-201-KD, BSM-201, BSM-202)
	Steel & Stainless Steel	(BSM-102-GS)
	ABS	(BMS-101-ABS)
Top :	Made of plywood	(BSM-102-GS, BMS-202)
	Made of Stainless Steel	(BSM-201-KD, BSM-201)
Food tray :	Molded plywood	(BSM-102-GS, BMS-202)
	Made of Stainless Steel with the bending process	(BSM-201-KD, BSM-201)
Finishing :	Epoxy Powder Coating & polished	(BSM-201-KD, BSM-201, BSM-202)
	Polished	(BSM-102-GS)
Mover :	By Gaspring	(BSM-102-GS)
	Manual adjuster	(BSM-201-KD, BSM-201, BMS - 101 - ABS, BMS-202)
Wheels :	2" Castor, two wheels wearing lockdown.	
Max Load :	15 Kg	
	size :	(L) 70 x (W) 40 x (H) 85 - 130 cm
Accessories :	Equipped with a locking pole lap (BSM-201-KD, BSM-201, BMS-101, BMS-202)	
Note :	Height and slope of the tray can be set with a round handle and lock (BSM-102-GS)	

Medicine Cabinet

BMC - 203 - 2P



BMC - 103 - 2P





BMC - 102 - 1P



BMC - 101 - 1P



BMC - 201 - 1P




Specification

Construction :	Steel Hollow, Plate is a material qualified and selected.	
Material :	Stainless Steel	(BMC-203-2P & BMC-201-1P)
	Steel	(BMC-103-2P, BMC-101-1P, BMC-102-1P)
Shelves :	3 Shelves Made of Glass	
Braket :	Made of Aluminium, retention bracket height can be adjusted according to need	
Type door :	1 Door	(BMC-201-1P, BMC-101-1P, BMC-102-1P)
	2 Door	(BMC-203-2P, BMC-103-2P)
Finishing :	Epoxy Powder Coating (BMC-103-2P, BMC-101-1P, BMC-102-1P) Polished (BMC-203-2P & BMC-201-1P)	
Side Guard :	Made of Glass	(BMC-103-2P, BMC-102-1P, BMC-203-2P, BMC-201-1P)
	Made of Plate	(BMC-101-1P)
Max Load :	20 Kg	
	size :	(L) 90 x (W) 45 x (H) 170 cm (BMC-203-2P, BMC-103-2P) (L) 70 x (W) 45 x (H) 170 cm (BMC-201-1P, BMC-101-1P, BMC-102-1P)
Accessories :	Lock door, Handle, Bracket Aluminium.	
Note :	The side of the cabinets are made of steel plate, shelves made of glass 5mm and adjustable, equipped with lock and handle. Under cabinets and 1 & 2 doors are made of steel plate fitted with a lock and handle	

Mortuary Refrigerator 2

BMO - 203 - 2PR

Specification




Construction :	Stainless Steel Hollow Pipe, Iron As, Plate, is a material qualified and selected.
Material :	Stainless steel hairline 304
Base :	Made of Stainless Steel
Finishing :	Hairline
Tray corpse :	Made of Stainless Steel
Wheels :	4", 4 -wheel Central locks
Max Load :	200 Kg
Size :	(L) 240 x (W) 90 x (H) 190 cm
Accessories :	Tray corpse, 1 Hp techumseh, refrigerant R 22
Power Consumption :	700 Watt
Power Supply :	220V - 50/60 Hz
Temperature :	0°C / -10°C
Note :	capacity for 2 bodies, there is a corpse tray, the temperature can be set to -10 °C, mobility is supported by the wheels.
Optional :	Strecher Cropse

Mortuary Carriage

BMO - 201-R

Specification




Construction :	Stainless Steel Pipe, Plate is a material qualified and selected.
Cover Mortuary :	Made of Stainless Steel Pipe & Cloth
Base :	Made of Stainless Steel Plate
Finishing :	Polished
Side Rail :	Made of Stainless Steel Pipe
I.V Stand :	-
Wheels :	4" Wheels
Max Load :	150 Kg
Weight :	- Kg
size :	(L) 200 x (W) 75 x (H) 80 cm
Accessories :	Cover Cloth
Note :	the cover using a hinge to facilitate opening and closing and locking features, for use under the front of the pipe and rubber serves as the locking movement and the rear wheels for ease of mobility.
Optional :	-

Mortuary Table

BMO - 201 - TR

Specification




Construction :	Stainless Steel Pipe, Hollow, Plate is a material qualified and selected.
Pillow the corpse :	Pillow foam semi latex covered with vinyl water repellent.
Finishing :	Polished
Side Rail :	-
I.V Stand :	-
Wheels :	-
Max Load :	150 Kg
Weight :	- Kg
size :	(L) 200 x (W) 70 x (H) 80 cm
Accessories :	Pillow the corpse
Note :	in order to facilitate and maintain the cleanliness of the corpse the table is equipped with cushions jsebagai bath mat, all stainless steel construction, there is sewerage.
Optional :	-

Operating Table Electric

BOT - 409 - 1ACS

Specification



Construction :	Stainless Steel Pipe, Hollow, Iron As, Plate, Steel Plate SPHC is a material qualified.
Base Plate :	Made of Steel Plate (SPHC 3mm)
Mattress :	Mattress foam semi latex covered with vinyl
Finishing :	Polished & Epoxy Powder Coating
Mover :	By Electrical Actuator Motor
Wheels :	4" & 3" wheel Central locks
Max Load :	200 Kg
Weight :	100 Kg
size :	(L) 204 x (W) 50 x (H) 85 - 105 cm
Accessories :	retaining shoulders, waist anchoring, anchoring hand, cloth hanger, retaining foot, Leg Support, Dirt Basin.
Position :	Tilting, High Low, Back Raise, Knee Raise, Trendeleunbrug & Anti Trendeleunbrug.
Power & Input Voltage :	200 Watt, 220V - 50/60 Hz
Note :	Actuator sistem 4000 - 6000 N, Gaspring 200 & 400N, Electrical use remote control, legs are divided into two parts.


Operating Table Electric

BOT - 210 - 2A/D

Specification



Construction :	Stainless Steel Pipe,Hollow, Iron As, Plate, Steel Plate SPHC is a material qualified.
Base Plate :	Made of Steel Plate (SPHC 3mm)
Mattress :	Mattress foam semi latex covered with vinyl
Finishing :	Polished & Epoxy Powder Coating
Mover	By Electrical Actuator Motor
Wheels	4" Suppo, 4 -wheel Central locks
Max Load :	200 Kg
Weight :	110 Kg
size :	(L) 204 x (W) 56 x (H) 85 - 105 cm
Accessories :	retaining shoulders, waist anchoring, anchoring hand, cloth hanger, retaining foot, Leg Support, Dirt basin.
Position :	Tilting, High Low, Back Raise,Knee Raise, Trendeleunbrug & Anti Trendeleunbrug.
Power & Input Voltage :	200 Watt, 220V - 50/60 Hz
Note:	Actuator sistem 4000 - 6000 N, 2 Units of energy storage batteries (can store electrical energy), Gastpring 200 & 400N, Electrical use remote control




Radiology Table

BRT - 401 - S

Specification



Construction :	Steel Hollow, Plate is a material qualified and selected.
Base Plate :	Made of Plywood & melaminto
Finishing :	Epoxy Powder Coating
Max Load :	150 Kg
Weight :	68,5 Kg
Wheel :	4" Castor, two wheels wearing lockdown.
size :	(L) 200 x (W) 70 x (H) 71 cm
Accessories :	film storage
Note :	there is a movie clip in storage, can be elongated handle is pulled out in a way that there is at storage of the film, in order to facilitate the movement of the film storage supported with sliding wheels. there is a rubber bottom feet for protection.
Optional :	-



Scrubbe Station Manual




BSR - 202 - 3 P



BSR - 402 - 2 PS



BSR - 401 - 1 PS



Specification

Construction :	Stainless Steel Hollow , pipe,Iron As,Plate is a material qualified and selected.
Body Cover :	made of Stainless Steel Plate
Finishing :	Polished
Light Poles :	Stainless Steel Plate.
Size:	(L) 90 x (W) 65 x (H) 180 cm (BSR-201-1PS)
	(L) 145 x (W) 65 x (H) 180 cm (BSR-202-2PS)
	(L) 260 x (W) 65 x (H) 150 cm (BSR-202-3P)
Power Consumption :	250 Watt
Power Supply :	AC 220 Volt, 50/60 Hz
Accessories :	Soap Tap, Water Tap, Hand Dryer & Mirror
system :	Manual
Sterilizing :	Filtering & U.V Sterilizer system
Filter :	Main filter 0.2 micron
Lamp :	TL U.V 30 watt Philips
Note:	using automated and automated manual system , there is a control panel boxes, door access
Optional :	The products can be adjusted to the wishes and needs of customers

Scrub Station Automatic




BSR - 202 - 3R



BSR - 402 - 2PR



BSR - 401 - 1PR



Specification

Construction :	Stainless Steel Hollow , pipe, Iron As, Plate is a material qualified and selected.	
Body Cover :	made of Stainless Steel Plate	
Finishing :	Polished	
Light Poles :	Stainless Steel Plate.	
Size:	(L) 90 x (W) 65 x (H) 180 cm	(BSR-401-1PR)
	(L) 145 x (W) 65 x (H) 180 cm	(BSR-402-2PR)
	(L) 260 x (W) 65 x (H) 150 cm	(BSR-202-3R),
Power Consumption :	250 Watt	
Power Supply :	AC 220 Volt, 50/60 Hz	
Accessories :	Soap Tap, Water Tap, Hand Dryer & Mirror Automatic & Manual	
Sterilizing :	Filtering & U.V Sterilizer system	
Filter :	Main filter 0.2 micron	
Lamp :	TL U.V 30 watt Philips	
Note:	using automated and automated manual system , there is a control panel boxes, door access	
Optional :	The products can be adjusted to the wishes and needs of customers	

Sound Proof Chamber



BSA - 401 - M

Specification

Construction :	Plywood, Glass, Steel plate, egg foam, gaspul is a material qualified and selected.
Base :	Made of Plywood Covered Gaspul
Finishing :	Melamic
Door :	made of plywood and foam egg gaspul
Handle :	Stainless Steel
Max Load :	100 Kg
Weight :	110 Kg
size :	(L) 90 x (W) 90 x (H) 200 cm
Volume Space :	(L) 80 x (W) 80 x (H) 175 cm
Accessories :	jack, door locks and handles, lamp, Shelves
Input Voltage :	220 Volt, 50/60 Hz
Note :	Frame made of plywood, foam absorbers made gaspul and eggs, has 5 control jack, 1 unit of light power button, input & out put,
System :	wearing knock down system.
Optimal :	-




Sound Proof Chamber



BSA - 401 - T

Specification

Construction :	Plywood, Glass, Steel plate, egg foam, gaspul Tacon is a material qualified and selected.
Base :	Made of Plywood Covered Gaspul
Finishing :	Tacon
Door :	Made of plywood and foam egg gaspul
Handle :	Stainless Steel
Max Load :	100 Kg
Weight :	120 Kg
size :	(L) 90 x (W) 90 x (H) 200 cm
Volume Space :	(L) 80 x (W) 80 x (H) 175 cm
Accessories :	jack, door locks and handles, lamp, Shelves
Input Voltage :	220 Volt, 50/60 Hz
Note :	Frame made of plywood, foam absorbers made gaspul and eggs, has 5 control jack, 1 unit of light power button, input & out put,
System :	wearing knock down system.
Optimal :	-



Stand Wascom



BSW - 201 - D



BSW - 101 - D



BSW - 201 - DD



BSW - 101 - S



BSW - 201 - S



BSW - 201 - SD

Specification


Construction :	Stainless Steel Pipe, Plate is a material qualified and selected.	
Material :	Steel	(BSW-101-D,BSW-101-S)
	Stainless Steel	(BSW-201-D,BSW-201-DD,BSW-201-S,BSW-201-D)
Base Plate :	Made of Stainless Steel plate	
Finishing :	Polished	(BSW-201-D, BSW-201-DD, BSW-201-S,BSW-201-SD)
	Epoxy Powder Coating	(BSW-101-D,BSW-101-S)
Type :	Single	(BSW-101-S, BSW-201-S, BSW-201-SD)
	Double	(BSW-201-D, BSW-201-DD)
Side Rail :	pipe	
Wheels :	2", two wheels with locking	
Max Load :	30 Kg	
Size :	(L) 50 x (W) 50 x (H) 80 cm	(BSW-101-S, BSW-201-S)
Note :	(L) 50 x (W) 40 x (H) 80 cm	(BSW-201-SD)
	(L) 75 x (W) 60 x (H) 80 cm	(BSW-101-D)
	(L) 75 x (W) 40 x (H) 80 cm	(BSW-201-D)
	(L) 110 x (W) 45 x (H) 80 cm	(BSW-201-DD)
Optional :	The detention basin washer is equipped with a towel. Stainless Steel plate (SPHC) as a solid resistance to cradle the basin.	
	Include basin	

Dressing Trolley



BDT - 201

BDT - 101



Specification


Construction :	Pipe,Hollow,Plate, iron is a material qualified and selected.	
Material :	Steel	(BDT-101)
	Stainless Steel	(BDT-201)
Base :	Made of Plate	
Finishing :	Epoxy Powder Coating	(BDT-101)
	Polished	(BDT-201)
Side Guard :	Iron Bar	
Handle :	Pipe	
Wheels :	3" , two wheels with locking	
Max Load :	50 Kg	
Size :	(L) 70 x (W) 40 x (H) 86 cm	
Accessories :	Handle, Drawer	
Note :	There are two shelves and a Iron as as side guards, equipped with handles and 2 storage drawer, there are 1 and bracket basin bucket.	
Optional :	Include basin & bucket	

Emergency Trolley S/S



BETR - 101

BETR - 201



Specification

Construction :	Stainless Steel Pipe, Hollow,Plate,Iron As, is a material qualified and selected.	
Material :	Steel	(BETR-101)
	Stainless Steel	(BETR-201)
Base :	Made of Plate	
Finishing :	Epoxy Powder Coating	(BETR-101)
	Polished	(BETR-201)
trays top :	Made of Stainless Steel Plate and hollow as well as side guard	
I.V Stand :	High-quality Stainless Steel material with 2 Infusion holder	
Wheels :	3" Wheel, two wheels with locking	
Max Load :	70 Kg	
size :	(L) 69 x (W) 50 x (H) 107 cm	
Accessories :	Infus Stand, Oxygen bracket,Electric connector Handle.	
Drawer :	there are 3 drawers made of plate	
Note :	materials are made of stainless steel & Steel, equipped with oxygen infusion stand locking knob star, 3-capacity electrical connector plug	
Optinal :	-	


Dirty Trolley

BFT - 208 - DFT

Specification



Construction :	Hollow , pipe, Iron Bar is a material qualified and selected
Material :	Stainless Steel
Finishing :	Polished
Under Prisoners :	Stainless Steel Iron As.
Body Frame :	High-quality Stainless Steel material
Wheels :	4 pcs 4" Wheel, 2 wheels wearing a safety lock
Dimensions :	(L) 91 x (W) 48 x (H) 110 cm
Accessories :	3 places dirty dishes and 2 trash
Note :	Frame is made of full steel, there are three places dirty dishes with plastic material and easy to clean, there are two trash for organic and non-organic.




USG Trolley Deluxe

Specification




BUT - 101 BUT - 101 - D

Construction :	Steel Hollow,Pipe, plate, Iron As, is a material qualified and selected.
Base :	Made of steel plate
Foundation tray :	-
Finishing :	Epoxy Powder Coating
Side Rail :	Made from a Iron as
I.V Stand :	-
Wheels :	3 " Wheel, two wheels with locking
Max Load :	50 Kg
Weight :	13 Kg
size :	(L) 58 x (W) 52 x (H) 101 cm
Accessories :	-
Note :	there are 3 shelves, top shelf lined with iron as in rool, the bottom shelf is made of a plate in bending and bounded as iron in the rool.
Optinal :	all materials are made of stainless steel




Food Trolley

Food Trolley Warmer 16 Tray
BFT - 202 - 16S



Food Trolley 16 Tray
BFT - 202 - 16S

Food Trolley 8 Tray
BFT - 201




Specification


Construction :	Stainless Steel Pipe, Hollow,Plate,Iron As, is a material qualified and selected.
Base :	Made of Stailless Steel Plate
Foundation tray :	Iron As in roll Stainless Steel Plate with bending
Finishing :	Polished
Side Rail :	Made of Stailless Steel Plate in bending
Wheels :	10" Wheel RHJTW, two wheels with locking 4" Wheel RHJTW, two wheels with locking
Max Load :	100 Kg
size :	(L) 90 x (W) 65 x (H) 120 cm
Accessories :	(L) 100 x (W) 55 x (H) 113 cm (L) 67 x (W) 57.5 x (H) 137 cm
Power Consumption :	Heater,Handle, Lock
Input Voltage :	600 Watt
Temperature :	220 Volt, 3A, 50/60 Hz
Door :	0 ^o - 110 °C
Note :	There are 1 & 2 doors made of Hollow and plate, which has two access doors. tray can accommodate as many as 16 racks, the top of the cart has a fence and a guard on the door there is a list of rubber to keep from going off the air.
Optinal :	Trays

Food Trolley


Food Trolley 24 Tray
BFT - 203



Food Trolley 32 Tray
BFT - 204 - 2F



Food Trolley 30 Tray
BFT - 211-30F




Specification


Construction :	Stainless Steel Pipe, Hollow, Plate, Iron As, is a material qualified and selected.	
Base :	Made of Stainless Steel Plate	
Foundation tray :	Stainless Steel Plate with bending	
Finishing :	Polished	
Side Rail :	Iron As in roll	
Wheels :	4" Wheel RHJTW, two wheels with locking	
Max Load :	100 Kg	
size :	(L) 100 x (W) 55 x (H) 160 cm	(BFT-203)
	L) 105 x (W) 55 x (H) 180 cm	(BFT-204-2F)
	L) 105 x (W) 65 x (H) 113 cm	(BFT-211-30F)
Accessories :	Handle, Lock door	
Door :	There are 4 doors made of Hollow and plate, which has two access doors front and rear	
Note :	tray can accommodate as many as 32 racks, trolley wearing the top of the fence guard and handle the driving	
Optinal :	Trays	

Instrument Trolley


Instrument Trolley P/C
BIT - 102




Instrument Trolley 2 Drawer
BIT - 102 - 1D



Instrument Trolley 2 Drawer
BIT - 101 - 2D



Instrument Glass Trolley P/C
BIT - 101



Specification

Construction :	Steel Pipe, Hollow, Plate is a material qualified and selected.	
Material :	Steel	
Base :	Made of Steel Plate	(BIT-102, BIT-102-1D, BIT-101-2D)
	Made of glass	(BIT-101)
Finishing :	Epoxy Powder Coating	
Side Guard :	Steel hollow	
Handle :	Steel Pipe	(BIT-102, BIT-102-1D, BIT-101-2D)
Wheels :	2" , two wheels with locking	
Max Load :	50 Kg	
Size :	(L) 60 x (W) 40 x (H) 80 cm	
Accessories :	Handle, Drawer	(BIT-102-1D, BIT-101-2D)
Note :	there are two shelves and a hollow pipe as side guards, equipped with handles	

Instrument Trolley



Clothes Trolley



Specification


Construction :	Pipe,Hollow,Plate is a material qualified and selected.	
Material :	Stainless steel	
Base :	Made of Stainless Steel Plate	(BIT-202, BIT-202-1D,BIT-201-2D,BIT-203)
	Made of glass	(BIT-201)
Finishing :	polished	
Side Guard :	Stainless Steel hollow	
Handle :	Stainless Steel Pipe (BIT-202, BIT-202-1D,BIT-201-2D,BIT-203)	
Wheels :	2" , two wheels with locking	
Max Load :	50 Kg	
Size :	(L) 70 x (W) 40 x (H) 80 cm	
Accessories :	Handle, Drawer	(BIT-202-1D,BIT-201-2D)
Note:	there are two shelves and a hollow pipe as side guards, equipped with handles	
	2 storage drawer.	(BIT-201-2D)
	1 storage drawer.	(BIT-202-1D)

Specification

Construction :	Stainless Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.	
Base :	Made of Steel sheet	(BLNT - 101, BLT - 101)
Finishing :	Epoxy Powder Coating	(BLNT - 101, BLT - 101)
	Polished	(BCT - 201) (BLNT - 201)
Side Rail :	iron as	(BLNT - 101)
	wire mesh	(BLT - 101)
Wheels :	4" wheel, 2 wheels wearing locks	
Max Load :	50 Kg	
	100 Kg	(BLNT - 101)
Size :	(L) 100 x (W) 55 x (H) 110 cm	(BLNT - 101)
	(L) 100 x (W) 55 x (H) 110 cm	(BLNT - 201)
	(L) 60 x (W) 60 x (H) 85 cm	(BCT - 201)
Accessories :	Linen	
Note :	part body made from steel plate in bending, the top has 2 shelves equipped with the side guide, the handle is made of stainless steel pipe. (BLNT - 101) (BLNT - 201)	
	frame is made of steel pipe are on a roll, then the cloth is made of strong linen and is impermeable to water. (BCT - 201)	
	there is a handle, the side of the trolley made of wiremesh (BLT - 101) (BLNT - 201)	
Optional :	all materials are made of stainless steel	

Clothes Trolley


Laundry Trolley
BLT - 210 - T



Laundr Cart
BCT - 202 - C



Laundry Trolley P/C
BLT - 101




Specification


Construction :	Stainless Steel Pipe,Hollow, Iron As, Plate is a material qualified and selected.	
Base :	Made of Steel sheet	(BLT - 210 - T, BLT - 101)
Finishing :	Epoxy Powder Coating	
	Polished	(BCT - 101) (BCT - 202 - C)
Side Rail :	iron as	(BCT - 202 - C)
	wire mesh	(BLT - 101)
Wheels :	4" wheel, 2 wheels wearing locks	
Max Load :	50 Kg	
	100 Kg	(BLT - 101)
Size :	(L) 76 x (W) 60 x (H) 89 cm	(BCT - 202 - C)
	(L) 105 x (W) 65 x (H) 80 cm	(BLT - 101)
	(L) 105 x (W) 65 x (H) 80 cm	(BLT - 210 - T)
Accessories :	Linen	
Note :	part body made from steel plate in bending, the top has 2 shelves equipped with the side guide, the handle is made of stainless steel pipe. (BLNT - 101)	
	frame is made of steel pipe are on a roll, then the cloth is made of strong linen and is impermeable to water. (BCT - 101)	
	there is a handle, the side of the trolley made of wiremesh (BLT - 101)	
Optional :	all materials are made of stainless steel	

Medicine Trolley


Medicine Trolley 24 Rak
BMT - 101 - 24R




Medicine Trolley 24 Rak
BMT - 103 - 24PL



Medicine Trolley 36 Rak
BMT - 101 - 36R



Medicine Trolley 36 Rak
BMT - 103 - 36PL



Specification


Construction :	Stainless Steel pipe, Steel Pipe,hollow,Plate, is a material qualified and selected.	
Material :	Steel	
Base :	Made of Steel Plate	
Finishing :	Epoxy Powder Coating	
Medicine box :	Made of Plastic	
	Made of Acrylic	
Handle :	Stainless Steel pipe with polished	
Wheels :	4" Wheel, two wheels with locking	
Max Load :	50 Kg	
	size :	(L) 51 x (W) 67 x (H) 100 cm (BMT-101-24R)
		(L) 70 x (W) 67 x (H) 100 cm (BMT-101-36R)
		(L) 80 x (W) 42 x (H) 104 cm (BMT-103-24PL)
		(L) 80 x (W) 42 x (H) 120 cm (BMT-103-36PL)
Accessories :	1 bottom drawers, door locks and handles	
	2 bottom drawers, door locks and handles	
Drawer :	Drawers made of plastic container	
Note :	material made of steel plate, bottom drawer is made of plate with the bending process.	
Optinal :	Material Stainless Steel	

Medicine Trolley


Medicine Trolley 30 Rak
BMT - 202 - 30R



Medicine Trolley 26 Rak
BMT - 201 - 26R



Medicine Trolley 24 Rak
BMT - 201 - 24R




Specification


Construction :	Stainless Steel Hollow,Plate,Stainless Steel Pipe, Acrylic is a material qualified and selected.	
Material :	Stainless Steel	
Base :	Made of Stainless Steel Plate	
Finishing :	Polished	
Medicine box :	Made of Acrylic	(BMT - 201 -24R & BMT - 202 -30R)
	Made of Stainless Steel	(BMT - 201 -26R)
Handle Trolley :	Stainless Steel pipe	
Wheels :	3" Wheel, two wheels with locking	
size :	(L) 80 x (W) 50 x (H) 80 cm	(BMT - 201 -26R)
	(L) 70 x (W) 67 x (H) 110 cm	(BMT - 202 -30R)
	(L) 51 x (W) 67 x (H) 100 cm	(BMT - 201 -24R)
Drawer :	24 drawers medicine, 2 bottom drawers and handles (BMT - 201 -26R)	
	30 drawers medicine, 2 bottom drawers door locks and handles (BMT - 202 -30R)	
	24 drawers medicine, 1 bottom drawers and handles (BMT - 201 -24R)	
Note :	Material made of Stainless steel plate, bottom drawer is made of Stainless Steel plate with the bending process, there is a storage rack on the back.	
Handle Drawers :	Aluminium	
Optinal :	Number of drawers can be adjusted to the needs	

Oxygen Trolley


BOXT-101-CMP




BOXT - 101 - 1.5



BOXT - 201 - 6



BOXT - 101 - 6



Specification


Construction :	Steel pipe, Plate, is a material qualified and selected.	
Material :	Stainless Steel	(BOXT-201-6)
	Steel	(BOXT-101-6,BOXT-101-CMP, BOXT-101-1.5)
Base :	Made of plate aeser	
Finishing :	Epoxy Powder Coating	(BOXT-101-6,BOXT-101-CMP, BOXT-101-1.5)
	Polished	(BOXT-201-6)
Side Rail :	Made of steel chain	
Wheels :	2" Wheel	(BOXT-101-CMP, BOXT-101-1.5)
	3" & 5" Wheel	(BOXT-101-6,BOXT-201-6)
size :	(L) 45 x (W) 108 x (H) 110 cm	(BOXT-101-6,BOXT-201-6)
	(L) 30 x (W) 25 x (H) 110 cm	(BOXT-101-CMP, BOXT-101-1.5)
Accessories :	1 unit of oxygen tubes 1m ³ , 1 set of regulator, 1 set of nasal canula	
Note :	1 oxygen tube fitted units, 1 set of regulator,1 set of nasal canula, there are rubber feet on the brake as	
Optional :	-	

U.V Mobile Sterilizer

Specification



Construction :	Stainless Steel Hollow , pipe,Iron As,Plate is a material qualified and selected.
Box :	made of Stainless Steel Plate
Finishing :	Polished
Light Poles :	Stainless Steel Iron As.
Handle :	High-quality Stainless Steel material
Wheels :	2 "using 2 wheel lock (BUW-201-4L) 3 "using 2 wheel lock (BUW-401-5LS)
Weight :	8.5 Kg
Size :	(L) 45 x (W) 45 x (H) 120 cm
Power Consumption :	126 Watt (BUW-401-5LS) 120 Watt (BUW-201-4L)
Power Supply :	AC 220 Volt, 50/60 Hz
Accessories :	Automatic Timer, Lamp, Bumper (BUW-401-5LS) Automatic Timer, Lamp (BUW-201-4L) Lamp 30W 4pcs & 6W 1 pcs (BUW-401-5LS)
Note :	Lamp 30 watt 4pcs (BUW-201-4L) Tilt adjustable lamp, the lamp is also equipped with a self-timer, the handle also functions as regulator of the tilt holder.




Visual Chart Light Box

Specification



BSN - 301 BSN - 302 - R

Construction :	Steel Plate, Acrylic, Aluminium List is a material qualified and selected.
Screen :	Acrylic
Finishing :	Epoxy Powder Coating
Body framework :	aluminium & Steel Plate in bending process
size :	(L) 32.5 x (W) 7.5 x (H) 92.5 cm
System :	Manual Control
Power Consumption :	24 Watt
Input voltage :	AC 220 - 230 Volt, 50/60 Hz
Accessories :	Taxing lamp 4 x 6 Watts, 12 volt, 5 A
Output :	visual chart light box using a system equipped with hinges and locks.
Note :	Equipped with remote control (BSN-302-R)




U.V Mobile Sterilizer

Specification



Construction :	Stainless Steel Hollow , pipe,Iron As,Plate is a material qualified and selected.
Box :	made of Stainless Steel Plate
Finishing :	Polished
Light Poles :	Stainless Steel Iron As.
Handle :	High-quality Stainless Steel material
Wheels :	2 "using 2 wheels plain and 2 wheel lock
Weight :	8.5 Kg
Size :	(L) 40 x (W) 40 x (H) 120 cm
Power Consumption :	126 Watt, 120 Watt.
Power Supply :	AC 220 Volt, 50/60 Hz
Accessories :	Automatic Timer, Lamp Lamp 30W 4pcs & 6W 1 pcs Lamp 30 watt 4pcs
Note :	Tilt adjustable lamp, the lamp is also equipped with a self-timer, the handle also functions as regulator of the tilt holder.
System Control :	Equipped with remote control




X - Ray Barrier

BXR - 301

Specification



Construction :	Steel Hollow, Plate elbow, Plate FB,plywood, melaminto is a material qualified and selected.
Frame :	Made of Steel plate elbow
Screen :	Made of plywood,plate PB, Glass PB & Melaminto
Finishing :	Epoxy Powder Coating
Pole :	Made of Steel hollow
Wheels :	4 " Wheel RHJTW, two wheels with locking
size :	(L) 120 x (W) 50 x (H) 180 cm
Note :	screen is made of plywood and melamine coated plate pb with adhesive glue EHA bon, to glass using glass pb with thickness of 8 mm and a size of 15 x 20 cm,
Optional :	all materials are made of stainless steel



X - Ray Viewer


BXR - 102 - LED



BXR - 101 - LED



BXR - 101 - S



BXR - 101 - D




Specification

Construction :	Steel Plate, Acrylic, Aluminium List is a material qualified and selected.
Screen :	Acrylic
Finishing :	Epoxy Powder Coating (BXR-101-S & BXR-101-D) Polished/Hair line (BXR-102-LED & BXR-101-LED)
Body framework :	Steel Plate in bending process Stainless Steel plate in bending process
size :	(L) 72.5 x (W) 4 x (H) 50 cm (BXR-101-LED) (L) 44 x (W) 4 x (H) 50 cm (BXR-102-LED) (L) 40 x (W) 9 x (H) 50 cm (BXR-101-S) (L) 70 x (W) 9 x (H) 50 cm (BXR-101-D)
Power Consumption :	DC 3 A, 4 Watt (BXR-101-LED) DC 3 A, 4 Watt (BXR-102-LED) 32 Watt (BXR-101-S) 64 Watt (BXR-101-D)
Input voltage :	AC 220 - 230 Volt, 50/60 Hz
Accessories :	Remot control, LED 1x4 Watts, adaptor charger
Output :	12 volt, 5 A, Electric Ballasts.
Note :	There is a list that serves as a locking aluminum film that will be read, Wearing the remote control

Standard Cassete

BXS - 101 - KD

Specification



Construction :	Steel Pipe, Plate, plywood is a material qualified and selected.
Foot step :	Made of Steel frame, plywood vinyl covered
Base :	Made of plywood vinyl covered
Finishing :	Epoxy Powder Coating & Polished
Pole cornerstone :	Made of Stainless Steel Pipe
Film storage :	Stainless Steel Plate in bending process
Max Load :	100 Kg
Weight :	14 Kg
size :	(L) 65 x (W) 50 x (H) 180 cm
Accessories :	Film storage & lock
Note :	the storage size of the film could be set high and low of the film and the object, all parts using knock-down system
Optional :	all materials are made of stainless steel

Product Population

ALKES. SIGOMPUL KAB. HUMBANG HASUNDUTAN
 ALKES. HUTAGALUNG KAB. HUMBANG HASUNDUTAN
 BALAIM KESEHATAN PARU, JABAR
 BLUD PKM. CEMPAKA PUTIH-JAKPUS
 BLUD PKM. KEC. SENEN - JAKPUS
 BLUD PKM. GROGOL PETAMBURAN
 BLUD PKM. KEC. TANJUNG PRIOK
 BLUD PKM. KEMBANGAN
 BLUD PUSKESMAS KEC. CIPAYUNG
 DINAS KESEHATAN KAB. ACEH TENGAH
 DINAS KESEHATAN KAB. BOJONEGORO
 DINAS KESEHATAN KAB. BREBES
 DINAS KESEHATAN KAB. ENDE
 DINAS KESEHATAN KAB. KERINCI
 DINAS KESEHATAN KOTA BANDUNG
 DINAS KESEHATAN KOTA PALEMBANG
 DINAS KESEHATAN KOTA PARIAMAN
 DINAS KESEHATAN PROV. BALI
 DINKES KAB. BANYUASIN
 DINKES KAB. BERAU
 DINKES KAB. BOGOR
 DINKES KAB. BOJONEGORO
 DINKES KAB. GARUT
 DINKES KAB. GUNUNG SITOLI
 DINKES KAB. INDRAMAYU
 DINKES KAB. JOMBANG
 DINKES KAB. KEDIRI
 DINKES KAB. OGAN KOMERING ULU
 DINKES KAB. PEKALONGAN
 DINKES KAB. PEMALANG
 DINKES KAB. PURWAKARTA
 DINKES KAB. SUKOHARJO
 DINKES KAB. SUMBA TENGAH
 DINKES KAB. TABALONG
 DINKES KAB. TABANAN
 DINKES KAB. TANGGAMUS
 DINKES KAB. TANJUNG JABUNG BARAT
 DINKES KAB. TULUNGAGUNG
 DINKES KOTA PALEMBANG
 DINKES KOTA SEMARANG
 DINKES KOTA SUNGAI PENUH
 DINKES PKM. KEDAUNG KOTA DEPOK
 DINKES PROV. BALI
 DINKES PROV. RIAU
 DINKES, KAB. MALANG
 DINKES, PENGENDALIAN PENDUDUK
 DINKES. KAB. ACEH TAMIANG
 DINKES. KAB. LAMPUNG SELATAN
 KANKES PELABUHAN KELAS II PADANG
 KEJAKSAAN AGUNG RI
 PEMDA DKI, RSUD PASAR MINGGU
 PEMDA KAB. PESISIR SELATAN
 PEMDA KAB. SERANG
 PEMDA KOTA BANDUNG
 PEMDA KOTA PALEMBANG

PKM BUGUL KIDUL KOTA PASURUAN
 PKM CAKUNG
 PKM CIKARANG
 PKM CIPAYUNG
 PKM DUREN SAWIT
 PKM I CILONGKOK KAB. BANYUMAS
 PKM KABUNAN KAB. PEMALANG
 PKM KAPUK
 PKM KEC. KEMBANGAN
 PKM KRAMAT JATI
 PKM MAJEGAN
 PKM MENTENG
 PKM PETARUKAN, KAB. PEMALANG
 PKM PULOGADUNG
 PKM RAWALO KAB. BANYUMAS
 PKM. BAKTIRAJA KAB. HUMBANG HASUNDUTAN
 PKM. BALONG KAB. PONOROGO
 PKM. BANDAR KEDUNG MULYO KAB. JOMBANG
 PKM. BANGETAYU PEMDA KOTA SEMARANG
 PKM. BANYUMAS KAB. BANYUMAS
 PKM. BENDOSARI KAB. SUKOHARJO
 PKM. BOJONG GEDE KAB. BOGOR
 PKM. BONAN DOLOK KAB. HUMBANG HASUNDUTAN
 PKM. BREBES
 PKM. BREBES KAB. BREBES
 PKM. CABANG BUNGIN
 PKM. CIBARUSAH
 PKM. CIGUDEG KAB. BOGOR
 PKM. CIKUYA KAB. TANGGERANG
 PKM. GUNUNGSARI KAB. BOJONEGORO
 PKM. HUTAGALUNG KAB. HUMBANG HASUNDUTAN
 PKM. HUTAPAUNG KAB. HUMBANG HASUNDUTAN
 PKM. I KEMRANJEN
 PKM. JEBED DINKES KAB. PEMALANG
 PKM. JOHAR BARU JAKARTA PUSAT
 PKM. KAB. GUNUNG SITOLI
 PKM. KABUNAN KAB. PEMALANG
 PKM. KARANG BAHAGIA
 PKM. KARANG HARJA
 PKM. KARANG MULYA
 PKM. KARANG SAMBUNG KAB. KEBUMEN
 PKM. KARANGAWEN I KAB. DEMAK
 PKM. KEBOAN KAB. JOMBANG
 PKM. KEC. KEBAYORAN BARU JAKSEL
 PKM. KEC. CIPAYUNG JAKTIM
 PKM. KLAREYAN KAB. PEMALANG
 PKM. LOJI KAB. MAJALENGKA
 PKM. MANISRENGGO
 PKM. MANISRENGGO KAB. KLATEN
 PKM. MATITI KAB. HUMBANG HASUNDUTAN
 PKM. MUNGO KAB. LIMA PULUH KOTA
 PKM. NANGGUNG KAB. BOGOR
 PKM. ONAN GANJANG KAB. HUMBANG HASUNDUTAN
 PKM. PAGUYANGAN KAB. BREBES
 PKM. PAKKAT KAB. HUMBANG HASUNDUTAN
 PKM. PAMANDEGAN KAB. LEBAK
 PKM. PANCORAN MAS 19 KOTA DEPOK

Product Population

PKM. PARANGINAN KAB. HUMBANG HASUNDUTAN
 PKM. PARLILITAN KAB. HUMBANG HASUNDUTAN
 PKM. PARUNGKUDA KAB. SUKABUMI
 PKM. PEBAYURAN
 PKM. PORIS PELAWAD KOTA TANGERANG
 PKM. PURWOKERTO BARAT KAB. BANYUMAS
 PKM. RAWALO KAB. BANYUMAS
 PKM. ROWOSARI KAB. PEMALANG
 PKM. SAITNIHUTA KAB. HUMBANG HASUNDUTAN
 PKM. SIGOMPUL KAB. HUMBANG HASUNDUTAN
 PKM. SUBAH KAB. BATANG
 PKM. SUKA TENANG
 PKM. SUMBERLAWANG KAB. SRAGEN
 PKM. TAJURHALANG KAB. BOGOR
 PKM. TANJUNG BREBES - BUMI INDAH
 PKM. TANJUNG TEJA
 PKM. TARABINTANG KAB. HUMBANG HASUNDUTAN
 PKM. TARUMAJAYA
 PKM. WEDI KAB. KLATEN
 PKM. WIRANG KAB. TABALONG
 PKM. WONOSALAM I KAB. DEMAK
 POLITEKNIK KESEHATAN DENPASAR
 POLITEKNIK KESEHATAN TANJUNG PINANG
 RS AL ISLAM BANDUNG
 RS ANNA MEDIKA
 RS ANTAM MEDIKA
 RS ARIYA MEDIKA
 RS ARJAWINANGUN CIREBON
 RS BUDI LESTARI
 RS BUNDA
 RS CIPUTRA
 RS CITAMA
 RS DINDA
 RS FATIMAH CILACAP
 RS GADING PLUIT
 RS HARAPAN KITA
 RS HARUM
 RS HERMINA BEKASI
 RS KANKER DHARMAIS
 RS KARTIKA HUSADA
 RS KELUARGA KITA
 RS KESDAM JAYA
 RS KOJA
 RS MARDHI WALUYO
 RS MASMITRA
 RS MAYANG MEDICAL CENTER
 RS MEDIMAS
 RS MEILIA
 RS MITRA KELUARGA BEKASI
 RS MITRA KEMAYORA
 RS MM INDRAMAYU
 RS MMC JAMBI
 RS MUHAMMADIYAH SRUWENG
 RS PERTAMINA BINA MEDIKA
 RS PERTAMINA SORONG
 RS POLRI KRAMAT JATI
 RS PURI ASIH

RS RUMKITAL
 RS SANGATTA OCCUPATION HEALTH CENTER
 RS SANTO YUSUP
 RS SEKAR KAMULIAN KUNINGAN
 RS SIAGA RAYA
 RS SINTANALA
 RS UIN
 RS UKI
 RS. KAB. SLEMAN
 RS. KUSTA SITANALA TANGERANG
 RS. PARU CISARUA
 RS. UMUM KEC. SAWAH BESAR JAKPUS
 RSAU Dr ERNASAWAN ANTARIKSA
 RSB SUN MARIE
 RSI PONDOK KOPI
 RSIA ASIH
 RSIA BUDI ASIH
 RSIA BUDI JAYA
 RSIA KARUNIA KASIH
 RSIA MUTIARA BUNDA
 RSIA RESTU KASIH
 RSIA SEKAR
 RSIA SELARAS
 RSIA SELASIH
 RSJ Dr. SOEHARTO HEEDJAN
 RSKD PROV. MALUKU
 RSU SERANG
 RSU SUMEDANG
 RSUD AL IKHSAN
 RSUD BANTEN
 RSUD BEKASI
 RSUD BREBES
 RSUD KAB. BUTON TENGAH
 RSUD KARAWANG
 RSUD KEC. CILINCING
 RSUD KOTA CILEGON
 RSUD MUARA DUA
 RSUD PROV. SULAWESI TENGGARA
 RSUK KALIDERES
 RSUK KOJA, JAKUT
 RSUK MAMPANG
 RSUKEC. TANGGERANG
 RSUP FATMAWATI
 RSUP HASAN SADIKIN BANDUNG
 RSUP PERSAHABATAN
 RUMKIT BAHYANGKARA PALANGKARAYA
 UNIVERSITAS GAJAH MADA
 UPTD PKM. MELUR KOTA PEKANBARU
 UPTD PKM. BUMIAYU KAB. BREBES
 UPTD PKM. KUTASARI KAB. PURBALINGGA
 UPTD PKM. PORIS PELAWAD TANGGERANG
 UPTD PKM. SIDANGWANGI KAB. MAJALENGKA
 UPTD PKM. TANJUNG KAB. BREBES
 WARINGIN KURUNG KAB. SERANG

kontak & informasi



PT. BUMI INDAH PUTRA

Shoow Room :

Graha Bumi Indah - lantai 2
Jln. Inspeksi Kalimantan Kav. Agraria Blok E No.10
Duren Sawit Jakarta Timur 13440
Telp.: 021 - 82731750 / 82731806
Fax.: 021 - 82731834

Office / Workshop :

Alex Bangun Cipendawa Baru
Rt/Rw 003/002 Bantar Gebang
Kota Bekasi - Jawa Barat
Telp.: 021 - 82731750 / 82731806
Fax.: 021 - 82731834



<https://maps.app.goo.gl/xedKKWpCsMAD2HqGA>



bipmed hospital equipments



www.bipmed.co.id

PT. BUMI INDAH PUTRA

Show Room :

Garaha Bumi Indah - Lantai 2
Jln. Inspeksi Kalimalang Kav. Agraria Blok E No. 10
Duren Sawit - Jakarta Timur 13440
Telp. : 021 - 82731750 / 8273106
Fax. : 021 - 82731834
Email : infobipmed@buniindah.com

Office / Workshop :

Alex Bangun Cipendawa Baru
Rt/Rw 003/002 Bantar Gebang
Kota Bekasi - Jawa Barat
Telp. : 021 - 82731750 / 82731806
Fax. : 021 - 82731834
Email : infobipmed@buniindah.com