



PT BUMI INDAH PUTRA

COVID-19 SPECIAL PRODUCT

Exclusive Distributor of:



Medical Waste Sterilizer




Workshop ✕

Alex Bangun Cipendawa Baru Rt. 003 Rw. 002
Kel. Bantar Gebang, Kec. Bantar Gebang
Kota Bekasi - Jawa Barat
Telp. 021 - 8273 1750 / 021-8273 1806
Fax. 021 - 8273 1834

Show Room 🏠

Graha Bumi Indah - Lantai 1
Jln. Inspeksi Kalimalang Kav. Agraria Blok E No. 10
Duren Sawit - Jakarta Timur 13430
Telp. 021 - 8660 8000 - ext 109, 114, 120
Fax. 021 - 8660 7111

**SWAB CHAMBER
BSC - 103 - SC**



**DISINFECTION CHAMBER
BDC - 102 - DC**



PRODUCT SPESIFICATION

Wall Frame	Steel hollow 30x30 mm
Material	Acrylic 5 mm transparent, Shark board 10 mm (PVC food grade)
Dimension	120 (L) x 120 (W) x 220 cm (H)
Power Consumption	70 Watt
Note	Speaker sound 2 way, light with TL, HEPA filter, & PVC gloves



Management System
ISO 9001:2008
www.tuv.com
ID 9105073007



PRODUCT SPESIFICATION

Construction	Steel hollow
Material	Acrylic 3 mm, steel hollow 30x30 plywood, & transparent curtain
Dimension	110 (L) x 110 (W) x 240 (H)
Tank Capacity	50 L - in (reception), 10 L - out (disposal)
Power Consumption	70 Watt
Note	Automatic sensor, light with TL, exhaust fan, nozzle spray




Management System
ISO 9001:2008
www.tuv.com
ID 9105073007



**ISOLATION STRETCHER
BAS-502-IS**



**ISOLATION CHAMBER
BAS - 500 - IS**



PRODUCT SPESIFICATION

Construction	Stainless Steel Pipe \varnothing 1"x1.2mm, \varnothing 5/8"x1.2mm, Plate t:1mm, t:3mm, Iron As, Aluminium Casting material qualified and selected.
Base	Made of Stainless Steel Plate
Mattress	Mattress foam semi latek covered with vynil, 60 mm thickness
Finishing	Polished
Dimension	(L)190x(W)60x(H)88 cm, (Bendable) (L)190 x(W)60 x(H)88 cm
Power supply	30 Watt-12Volt
Load capacity	Up to 150 Kg
System Filter	Inlet Prefilter and outlet HEPA filter for 0.3um particles
Accesories	Battery and charger



Management System ISO 9001:2015
www.tuv.com ID 9108627873



PRODUCT SPESIFICATION

Construction	Stainless Steel Pipe \varnothing 1"x1.2mm, \varnothing 5/8"x1.2mm, Plate t:1mm, t:3mm, Iron As, Aluminium Casting material qualified and selected.
Base	Made of Stainless Steel Plate
Mattress	Mattress foam semi latek covered with vynil, 60 mm thickness
Finishing	Polished
Dimension	(L)190x(W)60x(H)88 cm, (Bendable) (L)190 x(W)60 x(H)88 cm
Power supply	30 Watt-12Volt
Load capacity	Up to 150 Kg
System Filter	Inlet Prefilter and outlet HEPA filter for 0.3um particles
	Dilengkapi dengan Negative Pressure System
Accesories	Battery and charger




Management System ISO 9001:2015
www.tuv.com ID 9108627873




**ISOLATION STRETCHER DELUXE
BAS - 701 - IS - DELUXE**




BAS - 701- IS



Construction	Stainless stell pipe 1"x1,2mm, 5/8"x1,2mm, plate t:1mm t:3mm, iron as, aluminium casting, acrylic 5 mm, material qualified & selected
Base	Made of stainless steel plate
Finishing	Polishing
Dimension	(PxLxT) 213 x 66 x 137 cm
Dimension chamber	(PxLxT) 190 x 66 x 52 cm
HEPA filter	1 pcs, 0,3um, UV lamp : 10 watt
Power supply	AC 190-220V, 50?60Hz, Battery DC 12V 600mA
Component	Main chamber, HEPA filter, Negative presure system, 4 pcs glove portals, 2 pcs utility portals, battery & charger, chamber bag, stretcher



Management System
ISO 9001:2008
www.tuv.com
ID 9105073007



PRODUCT SPESIFICATION	
Contructions	Acrylic 2 mm transparent, Hollow 30x30x1.2 SS, 30x15x1.2 SS, Stainless Steel Hollow 20x20x1.2, plate t:1mm, pipe Ø1" x 1.2
Mattress	60 mm thickness
Finishing	Polished
Dimension	(L) 210 x (W) 63 x (H) 75 cm
Power Supply	450 watt
Load Capacity	up to 150 kg
System Filter	Inlet Prefilter and outlet HEPA filter for 0.3um particles
Accessories	Latex gloves, UV lamp, Fan & Battery charge



Management System ISO 9001:2015
www.tuv.com ID 9108627873





SERTIFIKAT
CPAKB
NO. FK.01.04/VI/953/2018



Manufactured by
PT. BUNDAH PUTRA

Portable Hepa Filter

BPH-102



High Quality Components and Filters :



SPESIFICATIONS

Dimension	: BPH-102: 51x32x89 cm
Air Volume	: 1000 - 1500 m3/h
Air Velocity	: 0.45-0.75 m/s
Material	: Stainless Steel
Power Consumption	: 60-125 w
Electric Standard	: 220 V AC, 50 Hz
Filtration Parts	: Prefilter, Hepa Filter, and UV sterilizer lamp
Noise	: 65 (H) db
HEPA Efficiency	: 99,97% 0.3 um
Replacement of filter	: Prefilter, Washable, per 3 months
HEPA Filter	: 6-12 months
Carbon Filter	: 6 months



Management System ISO 9001:2015
www.tuv.com ID 9108627873



Portable Hepa Filter

BPH-101



High Quality Components and Filters :



SPESIFICATIONS


Dimension	: BPH-102: 51x32x89 cm
Air Volume	: 1000 - 1500 m3/h
Air Velocity	: 0.45-0.75 m/s
Material	: Powder Steel Coating
Power Consumption	: 60-125 w
Electric Standard	: 220 V AC, 50 Hz
Filtration Parts	: Prefilter, Hepa Filter, and UV sterilizer lamp
Noise	: 65 (H) db
HEPA Efficiency	: 99,97% 0.3 um
Replacement of filter	: Prefilter, Washable, per 3 months
HEPA Filter	: 6-12 months
Carbon Filter	: 6 months




Management System ISO 9001:2015
www.tuv.com ID 9108627873



**WHEEL CHAIR
BKR-204-CVD**



**PETI JENAZAH COVID 19
BMO - 102 - C**



PRODUCT SPESIFICATION

Construction	Steel Pipe, Plate is a material qualified and selected.
Chamber	acrylic 8 mm
cushion	Synthetic leather is supple and strong
Finishing	Epoxy Powder Coating
Side Rail	pipe and placemat hand
Mover	Wheel
Wheels	Rear wheels 21" and front wheels 8"
Max Load	120 Kg
Weight	25 Kg
size	(L) 115 x (W) 70 x (H) 140 Cm
System Filter	Inlet Prefilter & outlet HEPA filter for 0.3um particles
Power Supply	450 watt
Accesories	UV lamp, fan and battery



Management System
ISO 9001:2008
www.tuv.com
ID 9105073007

PRODUCT SPESIFICATION

Material	Multyplex, kayu solid , kain & lapisan baha anti air
Finishing	Plamir, Cat putih
Dimension	P 200 x L 60 x T 50 Cm



Management System
ISO 9001:2015
www.tuv.com
ID 9105073007

